


# Transferrin Turbilatex®

The detection  
and quantification  
of **fecal transferrin**  
is a complementary tool  
in the **early diagnosis**  
of **bowel cancer**.



CerTest  
BIOTEC

**CerTest Biotec, S.L.**

Pol. Industrial Río Gállego II · Calle J, Nº1  
50840, San Mateo de Gállego, Zaragoza (Spain)

Tel. (+34) 976 520 354

Fax (+34) 976 106 268

certest@certest.es | turbilatex@certest.es

[www.certest.es](http://www.certest.es)



# Transferrin

Transferrin, like hemoglobin, is an iron-binding protein of the free iron in blood. Its presence in feces can be linked with bleeding of the upper or lower intestinal tract, **as it remains stable after its passage through the gastrointestinal tract. This makes it the perfect complementary analyte in the detection of occult digestive bleeding.**

This is the reason why CerTest has designed the **Transferrin Turbilatex®** turbidimetric assay, which detects and quantifies human transferrin in the same sample collection vial as **FOB Turbilatex®**.

**CerTest Transferrin Turbilatex®** is a non-invasive immunological test which can be used as a complement in the detection of FOBT. Together with **FOB Turbilatex®**, it is a sensitive specific tool for the detection of neoplasms in the colon.

## Our provided kits:



**FOB Turbilatex®**




**Calprotectin Turbilatex®**



**Transferrin Turbilatex®**



**Pancreatic elastase Turbilatex®**



**H. pylori Turbilatex®**

References	Description
TL-022TF100ED	<b>Transferrin Turbilatex</b> Combo kit 100 det.
TL-022TF200ED	<b>Transferrin Turbilatex</b> Combo kit 200 det.
MST-0019U	<b>Universal Turbilatex</b> Sample Collection Tube