


Transferrin Turbilatex®

The detection
and quantification
of **fecal transferrin**
is a complementary tool
in the **early diagnosis**
of **bowel cancer**.



CerTest
BIOTEC

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)

Tel. (+34) 976 520 354

Fax (+34) 976 106 268

certest@certest.es | turbilatex@certest.es

www.certest.es



Transferrin

Transferrin, like hemoglobin, is an iron-binding protein of the free iron in blood. Its presence in feces can be linked with bleeding of the upper or lower intestinal tract, **as it remains stable after its passage through the gastrointestinal tract. This makes it the perfect complementary analyte in the detection of occult digestive bleeding.**

This is the reason why CerTest has designed the **Transferrin Turbilatex®** turbidimetric assay, which detects and quantifies human transferrin in the same sample collection vial as **FOB Turbilatex®**.

CerTest Transferrin Turbilatex® is a non-invasive immunological test which can be used as a complement in the detection of FOBT. Together with **FOB Turbilatex®**, it is a sensitive specific tool for the detection of neoplasms in the colon.

Our provided kits:



FOB Turbilatex®



Calprotectin Turbilatex®



Transferrin Turbilatex®



Pancreatic elastase Turbilatex®



H. pylori Turbilatex®

References	Description
TL-022TF100ED	Transferrin Turbilatex Combo kit 100 det.
TL-022TF200ED	Transferrin Turbilatex Combo kit 200 det.
TL-022UN03E	Universal Turbilatex Sample diluent