



TECHNICAL SHEET

PRODUCT NAME: Human Elastase CELA3B recombinant protein

Code:	MT-25CELA3B, MT-25CELA3Bs1
Physical State:	Liquid
Source:	<i>Mammalian cells.</i>
Description:	Human Elastase CELA3B recombinant protein expressed and purified from <i>Mammalian cells</i> .
Appearance:	Clear
Purity:	>95%
Storage Conditions:	Short term storage at 2-8°C. Long term storage at -20°C. Avoid multiple/freeze thaw cycles by storing multiple aliquots at -20°C.
Health & Safety Information:	This product contains sodium azide as preservative. Although the amount is 0.09%, this substance is poisonous and hazardous and it should be handled by trained staff only. Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.
Presentation:	Tris buffer with additive and salts, pH 8.0; 0.09% Sodium Azide (NaN ₃).

THIS PRODUCT IS INTENDED FOR *IN VITRO* USE ONLY.

Date: 10/08/2018