



TECHNICAL SHEET

PRODUCT NAME: Human IL-6 recombinant protein

Code: MT-25IL6, MT-25IL6s1

Physical State: Liquid

Source: *Escherichia coli*

Description: Human IL-6 recombinant protein expressed and purified from *Escherichia coli*.

Appearance: Clear

Purity: >90%

Storage Conditions: Short term storage at 2-8°C. Long term storage at -20°C. Avoid multiple/freeze thaw cycles by storing multiple aliquots at -20°C.

Health & Safety Information: Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.

Presentation: Tris buffer with additive and salts, pH 7.5; 0.09% Sodium Azide (NaN₃).

THIS PRODUCT IS INTENDED FOR *IN VITRO* USE ONLY.

Date: 10/08/2018