



TECHNICAL SHEET

PRODUCT NAME: Anti-pancreatic Elastase mAb (clone EL04)

Code:	MT-16EL04
Clone:	EL04
Isotype:	IgA
Physical State:	Liquid
Host Animal (source):	Mouse
Specificity:	EL04 is specific for pancreatic Elastase detection.
Cross Reactivity:	No cross reactivity with pig elastase, bovine, pig and human haemoglobin, bovine, pig and human transferrin, human calprotectin, bovine and human lactoferrin.
Purification:	Anti-pancreatic Elastase mAb (clone EL04) has been purified in a multi-step procedure. ≥95% pure by SDS-PAGE.
Storage Conditions:	Storage at 2-8°C.
Health & Safety Information:	This product contains sodium azide as preservative. Although the amount is 0.09%, this substance is poisonous and hazardous and it should be handled by trained staff only. Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.
Presentation:	PBS pH: 7.4, 0.09% Sodium Azide (NaN ₃).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 23/10/2019