


Rapid Test

One-step immunochromatographic tests

Fast | Easy | Reliable | Simple



certest

IMMUNODIAGNOSTIC

Certest Immunodiagnostic is committed to providing a comprehensive portfolio of diagnostic reagents for detecting infectious diseases and different markers, covering different needs of formats, types of use, throughput, precision, and time from sample to result.


IMMUNOCHROMATOGRAPHY

provides rapid test solutions for aiding the professional of diagnostic with an immediate solution, needing no equipment, and having results usually in less than 15 minutes.

TURBIDIMETRY

provides diagnostics solutions for automated and quantitative results and higher throughputs.

CLIA provides a high sensitivity, automated diagnostic solution based on the chemiluminescent amplification of an immunocapture.



Rapid Test

Our products are the result of years of experience detecting the presence of antigens with the lateral flow **immunochromatographic** technique.

The devices are based on the coloured latex particles. Showing results on bands of different colours facilitates the interpretation of the result.

Certest Lateral Flow Immunochromatography portfolio mainly covers the diagnostic of the gastrointestinal, respiratory infectious diseases and the detection of inflammatory or tumor markers in the Gastrointestinal tract.

Certest advantages


- ✔ **Single device**, multiple tests from one sample.
- ✔ **One sample**, one preparation, panel testing.
- ✔ Multiple analytes diagnosed in **10 minutes**.
- ✔ **Reduced risk** of sample and result mishandling.
- ✔ Operational efficiency.



Multianalyte COMBO
(1, 2, 3 or 4 strips)

Our cut-off tests are equipped with special vials which facilitates fecal sample collection.

Results in
10 minutes



Certest ResisCheck® Carbapenemases

- ✔ **IVDR Certification**.
- ✔ **European manufacturing**.
- ✔ **High signal intensity** across all 5 lines.
- ✔ **Clear results** with a clean membrane background.
- ✔ **Larger detection window** for reading results.
- ✔ **Differential control line** in green.
- ✔ **Sample type:** Bacterial isolates from any biological specimen.
- ✔ **Results in 15 minutes**.



Fecal antigens	Respiratory antigens
 <p>VIRUSES</p> <ul style="list-style-type: none"> Adenovirus  Astrovirus  Enterovirus Norovirus GI/GII  Rotavirus  <p>BACTERIA</p> <ul style="list-style-type: none"> Helicobacter pylori  Clostridium difficile antigen GDH  Escherichia coli O157 Campylobacter  Salmonella Salmonella Typhi Salmonella Paratyphi Yersinia O:3 Yersinia O:9 <p>PARASITES</p> <ul style="list-style-type: none"> Cryptosporidium parvum  Giardia lamblia  Entamoeba  	<p>VIRUSES</p> <ul style="list-style-type: none"> Flu A RSV Adenovirus Resp. SARS-CoV-2 <p>BACTERIA</p> <ul style="list-style-type: none"> Strep. A
 <p>VIRUSES</p> <ul style="list-style-type: none"> Rotavirus + Adenovirus  Norovirus GI + Norovirus GII  Rotavirus + Norovirus  <p>BACTERIA</p> <ul style="list-style-type: none"> H. pylori + Transferrin Clostridium difficile (GDH) + Lactoferrin Clostridium difficile Toxin A + B  EHEC VT1 + VT2 Salmonella typhi + paratyphi Yersinia enterocolitica O:3 + O:9 <p>PARASITES</p> <ul style="list-style-type: none"> Crypto + Giardia  	<p>VIRUSES</p> <ul style="list-style-type: none"> Flu A + Flu B RSV + Adenovirus Resp.
 <p>VIRUSES</p> <ul style="list-style-type: none"> Rotavirus + Adenovirus + Astrovirus  Rotavirus + Adenovirus + Norovirus  <p>BACTERIA</p> <ul style="list-style-type: none"> Clostridium difficile GDH + Toxin A + Toxin B  E. coli O157 + VT1 + VT2 <p>PARASITES</p> <ul style="list-style-type: none"> Crypto + Giardia + Entamoeba  	<p>VIRUSES</p> <ul style="list-style-type: none"> Flu A + Flu B + RSV SARS-CoV-2 + Flu A + Flu B
 <p>VIRUSES</p> <ul style="list-style-type: none"> Rotavirus + Adenovirus + Astrovirus + Norovirus  <p>BACTERIA</p> <ul style="list-style-type: none"> C. Difficile GDH + Toxin A + Toxin B + Lactoferrin 	<p>VIRUSES</p> <ul style="list-style-type: none"> Flu A + Flu B + RSV + Adenovirus Resp. SARS-CoV-2 + Flu A + Flu B + RSV

Inflammatory and tumor markers



MARKERS

Hemoglobin (FOB)
Transferrin
Calprotectin
Lactoferrin

Urine

BACTERIA

Legionella
Streptococcus pneumoniae

MARKERS



FOB 50 + 200
FOB + Transferrin
H. pylori + Transferrin
Calprotectin 50 + 200
Calprotectin + Lactoferrin
Clostridium difficile (GDH) + Lactoferrin

BACTERIA

Streptococcus pneumoniae + Legionella

MARKERS



FOB + Transferrin + Calprotectin + Lactoferrin
GDH + Toxin A + Toxin B + Lactoferrin

Antimicrobial resistance



KPC | NDM | VIM | OXA-48 | IMP



Our tests have two differentiated lines.

Control line - green color
Positive line - red color



Possibility to add
positive and negative controls.





 **Fecal Antigens**

 **GASTROINTESTINAL VIRUS**

 **ADENOVIRUS**

Reference	Product name
AD820001V	Adenovirus (20 test cassettes with stool vials) 
AD820001VC	Adenovirus (20 test cassettes with stool vials + Positive Control) 
AD220001P	Adenovirus 20 Positive Controls Kit



 **ASTROVIRUS**

Reference	Product name
AT820001V	Astrovirus (20 test cassettes with stool vials) 
AT820001VC	Astrovirus (20 test cassettes with stool vials + Positive Control) 
AT220001P	Astrovirus 20 Positive Controls Kit



 **ENTEROVIRUS**

Reference	Product name
EV820001V	Enterovirus (20 test cassettes with stool vials)
EV820001VC	Enterovirus (20 test cassettes with stool vials + Positive Control)
EV220001P	Enterovirus 20 Positive Controls Kit







 **NOROVIRUS**

Reference	Product name
NN820001V	Norovirus (20 test cassettes with stool vials) 
NN820001VC	Norovirus (20 test cassettes with stool vials + Positive Control) 
NP220001P	Norovirus 20 Positive Controls Kit





 **ROTAVIRUS**

Reference	Product name	
R820001V	Rotavirus (20 test cassettes with stool vials)	
R820001VC	Rotavirus (20 test cassettes with stool vials + Positive Control)	
R220001P	Rotavirus 20 Positive Controls Kit	



 **COMBO**

Reference	Product name	
XR882001V	Rotavirus + Adenovirus (20 test combo cassettes with stool vials)	
XR882001VC	Rotavirus + Adenovirus (20 test combo cassettes with stool vials + Positive Controls)	
NV882001V	Norovirus GI + GII (20 test combo cassettes with stool vials)	
NV882001VC	Norovirus GI + GII (20 test combo cassettes with stool vials + Positive Controls)	
RN882001V	Rotavirus + Noro GI/GII (20 test combo cassettes with stool vials)	
RN882001VC	Rotavirus + Noro GI/GII (20 test combo cassettes with stool vials + Positive Controls)	

 **TRIPLE**

Reference	Product name	
RT872001V	Rotavirus + Adenovirus + Astrovirus (20 test triple cassettes with stool vials)	
RT872001VC	Rotavirus + Adenovirus + Astrovirus (20 test triple cassettes with stool vials + Positive Controls)	
RA872001V	Rotavirus + Adenovirus + Norovirus (20 test triple cassettes with stool vials)	
RA872001VC	Rotavirus + Adenovirus + Norovirus (20 test triple cassettes with stool vials + Positive Controls)	

 **QUAD**



Reference	Product name	
RG862001V	Rota + Adeno + Astro + Noro (20 test quad cassettes with stool vials)	
RG862001VC	Rota + Adeno + Astro + Noro (20 test quad cassettes with stool vials + Positive Controls)	





Fecal Antigens

GASTROINTESTINAL BACTERIA

HELICOBACTER PYLORI

Reference	Product name	
P820001V	H. Pylori (20 test cassettes with stool vials)	
P820001VC	H. Pylori (20 test cassettes with stool vials + Positive Control)	
P220001P	H. Pylori 20 Positive Controls Kit	



CLOSTRIDIUM DIFFICILE

Reference	Product name	
GD820001V	Clostridium Difficile antigen GDH (20 test cassettes with stool vials)	
GD820001VC	Clostridium Difficile antigen GDH (20 test cassettes with stool vials + Positive Control)	
GD220001P	Clostridium Difficile antigen GDH 20 Positive Controls Kit	
TA220001P	Clostridium Difficile Toxin A 20 Positive Controls Kit	
TB220001P	Clostridium Difficile Toxin B 20 Positive Controls Kit	

ESCHERICHIA COLI

Reference	Product name	
E820001V	E. Coli O157:H7 (20 test cassettes with stool vials)	
E820001VC	E. Coli O157:H7 (20 test cassettes with stool vials + Positive Control)	
E220001P	E. Coli O157:H7 20 Positive Controls Kit	
VT220001P	EHEC VT1 Verotoxin 1, 20 Positive Controls Kit	
VX220001P	EHEC VT2 Verotoxin 2, 20 Positive Controls Kit	

CAMPYLOBACTER

Reference	Product name	
CA820001V	Campylobacter (20 test cassettes with stool vials)	
CA820001VC	Campylobacter (20 test cassettes with stool vials + Positive Control)	
CA220001P	Campylobacter 20 Positive Controls Kit	

SALMONELLA

Reference	Product name
SA820001V	Salmonella (20 test cassettes with stool vials)
SA820001VC	Salmonella (20 test cassettes with stool vials + Positive Control)
SA220001P	Salmonella 20 Positive Controls Kit

SALMONELLA TYPHI

Reference	Product name
ST820001V	Salmonella Typhi (20 test cassettes with stool vials)
ST820001VC	Salmonella Typhi (20 test cassettes with stool vials + Positive Control)
ST220001P	Salmonella Typhi 20 Positive Controls Kit

SALMONELLA PARATYPHI

Reference	Product name
SP820001V	Salmonella Paratyphi (20 test cassettes with stool vials)
SP820001VC	Salmonella Paratyphi (20 test cassettes with stool vials + Positive Control)
SP220001P	Salmonella Paratyphi 20 Positive Controls Kit

YERSINISA O:3

Reference	Product name
YE820001V	Yersinia Enterocolitica O:3 (20 test cassettes with stool vials)
YE820001VC	Yersinia Enterocolitica O:3 (20 test cassettes with stool vials + Positive Control)
YE220001P	Yersinia Enterocolitica O:3 20 Positive Controls Kit

YERSINISA O:9



Reference	Product name
YC820001V	Yersinia Enterocolitica O:9 (20 test cassettes with stool vials)
YC820001VC	Yersinia Enterocolitica O:9 (20 test cassettes with stool vials + Positive Control)
YC220001P	Yersinia Enterocolitica O:9 20 Positive Controls Kit



Fecal Antigens



COMBO

Reference	Product name	
CD882001V	Clostridium Difficile Toxin A + Toxin B (20 test combo cassettes with stool vials)	
CD882001VC	Clostridium Difficile Toxin A + Toxin B (20 test combo cassettes with stool vials + Positive Controls)	
PT882001VF	H. Pylori + Transferrin (20 test combo cassettes with stool vials)	
PT882001VFC	H. Pylori + Transferrin (20 test combo cassettes with stool vials + Positive Controls)	
GF882001VF	GDH + Lactoferrin (20 test combo cassettes with stool vials)	
GF882001VFC	GDH + Lactoferrin (20 test combo cassettes with stool vials + Positive Controls)	
SX882001V	EHEC VT1 + VT2 (20 test combo cassettes with stool vials)	
SX882001VC	EHEC VT1 + VT2 (20 test combo cassettes with stool vials + Positive Controls)	
TP882001V	Salmonella Typhi + Paratyphi (20 test combo cassettes with stool vials)	
TP882001VC	Salmonella Typhi + Paratyphi (20 test combo cassettes with stool vials + Positive Controls)	
YS882001V	Yersinia Enterocolitica O:3 + O:9 (20 test combo cassettes with stool vials)	
YS882001VC	Yersinia Enterocolitica O:3 + O:9 (20 test combo cassettes with stool vials + Positive Controls)	



TRIPLE

Reference	Product name	
GX872001V	C. Difficile GDH + Toxin A + Toxin B (20 test triple cassettes with stool vials)	
GX872001VC	C. Difficile GDH + Toxin A + Toxin B (20 test triple cassettes with stool vials + Positive Controls)	
EX872001V	E. Coli O157:H7 + VT1 + VT2 (20 test triple cassettes with stool vials)	
EX872001VC	E. Coli O157:H7 + VT1 + VT2 (20 test triple cassettes with stool vials + Positive Controls)	





QUAD



Reference	Product name	
LX862001VF	C. Difficile GDH + Toxin A + Toxin B + Lactoferrin (20 test quad cassettes with stool vials)	
LX862001VFC	C. Difficile GDH + Toxin A + Toxin B + Lactoferrin (20 test quad cassettes with stool vials + Positive Controls)	

GASTROINTESTINAL PARASITES



CRYPROSPORIDIUM

Reference	Product name	
K820001V	Crypto (20 test cassettes with stool vials)	
K820001VC	Crypto (20 test cassettes with stool vials + Positive Control)	
K220001P	Crypto 20 Positive Controls Kit	



GIARDIA

Reference	Product name	
GW820001V	Giardia+ (20 test cassettes with stool vials)	
GW820001VC	Giardia+ (20 test cassettes with stool vials + Positive Control)	
GW220001P	Giardia+ 20 Positive Controls Kit	



ENTAMOEBIA

Reference	Product name	
EH820001V	Entamoeba (20 test cassettes with stool vials)	
EH820001VC	Entamoeba (20 test cassettes with stool vials + Positive Control)	
EH220001P	Entamoeba 20 Positive Controls Kit	

COMBO

Reference	Product name	
NW882001V	Crypto + Giardia (20 test combo cassettes with stool vials)	
NW882001VC	Crypto + Giardia (20 test combo cassettes with stool vials + Positive Controls)	

TRIPLE

Reference	Product name	
NE872001V	Crypto + Giardia + Entamoeba (20 test combo cassettes with stool vials)	
NE872001VC	Crypto + Giardia + Entamoeba (20 test combo cassettes with stool vials + Positive Controls)	



Inflammatory & Tumor Makers



FOB

Reference	Product name
F820001F	FOB (20 test cassettes with stool vials)
F820001FC	FOB (20 test cassettes with stool vials + Positive Control)
F220001P	FOB 20 Positive Controls Kit



CALPROTECTIN

Reference	Product name
CP820001F	Calprotectin (20 test cassettes with stool vials)
CP820001FC	Calprotectin (20 test cassettes with stool vials + Positive Control)
CP220001P	Calprotectin 20 Positive Controls Kit



TRANSFERRIN

Reference	Product name
TF820001F	Transferrin (20 test cassettes with stool vials)
TF820001FC	Transferrin (20 test cassettes with stool vials + Positive Control)
TF220001P	Transferrin 20 Positive Controls Kit



LACTOFERRIN

Reference	Product name
LT820001F	Lactoferrin (20 test cassettes with stool vials)
LT820001FC	Lactoferrin (20 test cassettes with stool vials + Positive Control)
LT220001P	Lactoferrin 20 Positive Controls Kit



COMBO

Reference	Product name
FB882001F	FOB 50 + 200 (20 test combo cassettes with stool vials)
FB882001FC	FOB 50 + 200 (20 test combo cassettes with stool vials + Positive Controls)

CC882001F	Calprotectin 50 + 200 (20 test combo cassettes with stool vials)
CC882001FC	Calprotectin 50 + 200 (20 test combo cassettes with stool vials + Positive Controls)
FT882001F	FOB + Transferrin (20 test combo cassettes with stool vials)
FT882001FC	FOB + Transferrin (20 test combo cassettes with stool vials + Positive Controls)
CL882001F	Calprotectin + Lactoferrin (20 test combo cassettes with stool vials)
CL882001FC	Calprotectin + Lactoferrin (20 test combo cassettes with stool vials + Positive Controls)



QUAD

Reference	Product name
FC862001F	FOB + Transferrin + Calprotectin + Lactoferrin (20 test combo cassettes with stool vials)
FC862001FC	FOB + Transferrin + Calprotectin + Lactoferrin (20 test combo cassettes with stool vials + Positive Controls)



Urine Antigens

URINE BACTERIA



LEGIONELLA PNEUMOPHILA

Reference	Product name
LP820201PC	Legionella (20 test cassettes with diluent bottle + Positive Control + Negative Control)
LP220001P	Legionella 20 Positive Controls Kit



STREPTOCOCCUS PNEUMONIAE

Reference	Product name
SN820201PC	S. pneumoniae (20 test cassettes with diluent bottle + Positive Control + Negative Control)
SN220001P	S. pneumoniae 20 Positive Controls Kit



COMBO

Reference	Product name
NL882201PC	S. pneumoniae + Legionella (20 test cassettes with diluent bottle + Positive Controls + N. C.)



Respiratory Antigens

RESPIRATORY VIRUS

INFLUENZA

Reference	Product name
Y820001PC	Influenza A (20 test cassettes with sample vials + Positive Control)
Y220001P	Influenza A 20 Positive Controls Kit
YB220001P	Influenza B 20 Positive Controls Kit

RSV (RESPIRATORY SYNCYTIAL VIRUS)

Reference	Product name
RV820001PC	RSV (20 test cassettes with sample vials + Positive Control)
RV220001P	RSV 20 Positive Controls Kit

ADENOVIRUS RESPIRATORY

Reference	Product name
AP820001PC	Adenovirus Resp. (20 test cassettes with sample vials + Positive Control)
AP220001P	Adenovirus Resp. 20 Positive Controls Kit

SARS COV-2

Reference	Product name
SC820001P	SARS CoV-2 (20 test cassettes with sample vials)
SC820001PC	SARS CoV-2 (20 test cassettes with sample vials + Positive Control)
SC860001PC	SARS CoV-2 (60 test cassettes with sample vials + Positive Control)
SC220001P	SARS CoV-2 20 Positive Controls Kit

COMBO

Reference	Product name
Z882001PC	Influenza A + B (20 test cassettes with sample vials + Positive Control)
VR882001PC	RSV + Adenovirus Resp. (20 test cassettes with sample vials + Positive Control)



TRIPLE

Reference	Product name
ZS872001PC	Influenza A + B + RSV (20 test cassettes with sample vials + Positive Control)
SS872001PC	SARS CoV-2 + Flu A + Flu B (20 test cassettes with sample vials + Positive Control)



QUAD

Reference	Product name
ZV862001PC	Influenza A + B + RSV + Adeno Resp. (20 test cassettes with sample vials + Positive Control)
SV862001PC	SARS CoV-2 + Flu A + Flu B + RSV (20 test cassettes with sample vials + Positive Control)

RESPIRATORY BACTERIA

STREP. A

Reference	Product name
T820201PC	Strep A (20 test cassettes with diluent bottles + tubes and pipettes vials + Positive Control)
T220001P	Strep A 20 Positive Controls Kit



Antimicrobial resistance

Reference	Product name
RB820001VC	Certest ResisCheck® Carbapenemases (20 test cassettes with sample vials + Liophilised Positive Control)



Any test KIT which includes positive control could have negative control too adding C at the end of the reference.

Reference	Product name
CN320001P	Negative Control 20 Negative Controls Kit

Certest Immunodiagnostic offers simple, sensitive, specific screening assays for a presumptive diagnosis of infectious diseases.

Advantages:



Noninvasive diagnostic



All included. No additional equipment needed



Low operational costs



Immediate results



Ease of use and interpretation

From Discovery
to Delivery,
Biotech you
can trust.

certest

IMMUNODIAGNOSTIC

Certest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, N°1
50840, San Mateo de Gállego, Zaragoza (Spain)
Tel. (+34) 976 520 354
sales@certest.es
www.certest.es



Find us on

@certestbiotec

Certest Rapid Test/GEN-0426EN

Modification rights reserved. All rights reserved. © Certest Biotec, S.L.
The products, services and data set out in this document may suffer changes and/or variations on the texts and pictures shown.