


# Rapid Test

One-step immunochromatographic tests

Fast | Easy | Reliable | Simple



**certest**

IMMUNODIAGNOSTIC

Certest Immunodiagnostic is committed to providing a comprehensive portfolio of diagnostic reagents for detecting infectious diseases and different markers, covering different needs of formats, types of use, throughput, precision, and time from sample to result.


### IMMUNOCHROMATOGRAPHY

provides rapid test solutions for aiding the professional of diagnostic with an immediate solution, needing no equipment, and having results usually in less than 15 minutes.

### TURBIDIMETRY

provides diagnostics solutions for automated and quantitative results and higher throughputs.

**CLIA** provides a high sensitivity, automated diagnostic solution based on the chemiluminescent amplification of an immunocapture.



## Rapid Test

Our products are the result of years of experience detecting the presence of antigens with the lateral flow **immunochromatographic** technique.

The devices are based on the coloured latex particles. Showing results on bands of different colours facilitates the interpretation of the result.


**Certest Lateral Flow Immunochromatography** portfolio mainly covers the diagnostic of the gastrointestinal, respiratory infectious diseases and the detection of inflammatory or tumor markers in the Gastrointestinal tract.

- ✔ **Single device**, multiple tests from one sample.
- ✔ **One sample**, one preparation, panel testing.
- ✔ Multiple analytes diagnosed in **10 minutes**.
- ✔ **Reduced risk** of sample and result mishandling.
- ✔ Operational efficiency.



**Multianalyte COMBO (1, 2, 3 or 4 strips)**

*Our cut-off tests are equipped with special vials which facilitates fecal sample collection.*



Fecal antigens	Respiratory antigens
 <p><b>VIRUSES</b></p> <ul style="list-style-type: none"> <li>Adenovirus </li> <li>Astrovirus </li> <li>Enterovirus</li> <li>Norovirus GI/GII </li> <li>Rotavirus </li> </ul> <p><b>BACTERIA</b></p> <ul style="list-style-type: none"> <li>Helicobacter pylori </li> <li>Clostridium difficile antigen GDH </li> <li>Escherichia coli O157</li> <li>Campylobacter </li> <li>Salmonella</li> <li>Salmonella Typhi</li> <li>Salmonella Paratyphi</li> <li>Yersinia O:3</li> <li>Yersinia O:9</li> </ul> <p><b>PARASITES</b></p> <ul style="list-style-type: none"> <li>Cryptosporidium parvum </li> <li>Giardia lamblia </li> <li>Entamoeba </li> </ul>	<p><b>VIRUSES</b></p> <ul style="list-style-type: none"> <li>Flu A</li> <li>RSV</li> <li>Adenovirus Resp.</li> <li>SARS-CoV-2</li> </ul> <p><b>BACTERIA</b></p> <ul style="list-style-type: none"> <li>Strep. A</li> </ul>
 <p><b>VIRUSES</b></p> <ul style="list-style-type: none"> <li>Rotavirus + Adenovirus </li> <li>Norovirus GI + Norovirus GII </li> <li>Rotavirus + Norovirus </li> </ul> <p><b>BACTERIA</b></p> <ul style="list-style-type: none"> <li>H. pylori + Transferrin</li> <li>Clostridium difficile (GDH) + Lactoferrin</li> <li>Clostridium difficile Toxin A + B </li> <li>EHEC VT1 + VT2</li> <li>Salmonella typhi + paratyphi</li> <li>Yersinia enterocolitica O:3 + O:9</li> </ul> <p><b>PARASITES</b></p> <ul style="list-style-type: none"> <li>Crypto + Giardia </li> </ul>	<p><b>VIRUSES</b></p> <ul style="list-style-type: none"> <li>Flu A + Flu B</li> <li>RSV + Adenovirus Resp.</li> </ul>
 <p><b>VIRUSES</b></p> <ul style="list-style-type: none"> <li>Rotavirus + Adenovirus + Astrovirus </li> <li>Rotavirus + Adenovirus + Norovirus </li> </ul> <p><b>BACTERIA</b></p> <ul style="list-style-type: none"> <li>Clostridium difficile GDH + Toxin A + Toxin B </li> <li>E. coli O157 + VT1 + VT2</li> </ul> <p><b>PARASITES</b></p> <ul style="list-style-type: none"> <li>Crypto + Giardia + Entamoeba </li> </ul>	<p><b>VIRUSES</b></p> <ul style="list-style-type: none"> <li>Flu A + Flu B + RSV</li> <li>SARS-CoV-2 + Flu A + Flu B</li> </ul>
 <p><b>VIRUSES</b></p> <ul style="list-style-type: none"> <li>Rotavirus + Adenovirus + Astrovirus + Norovirus </li> </ul> <p><b>BACTERIA</b></p> <ul style="list-style-type: none"> <li>C. Difficile GDH + Toxin A + Toxin B + Lactoferrin</li> </ul>	<p><b>VIRUSES</b></p> <ul style="list-style-type: none"> <li>Flu A + Flu B + RSV + Adenovirus Resp.</li> <li>SARS-CoV-2 + Flu A + Flu B + RSV</li> </ul>

**Inflammatory and tumor markers**

**Urine**



**MARKERS**

- Hemoglobin (FOB)
- Transferrin
- Calprotectin
- Lactoferrin

**BACTERIA**

- Legionella
- Streptococcus pneumoniae



**MARKERS**

- FOB 50 + 200
- FOB + Transferrin
- H. pylori + Transferrin
- Calprotectin 50 + 200
- Calprotectin + Lactoferrin
- Clostridium difficile (GDH) + Lactoferrin

**BACTERIA**

- Streptococcus pneumoniae + Legionella



**MARKERS**

- FOB + Transferrin + Calprotectin + Lactoferrin
- GDH + Toxin A + Toxin B + Lactoferrin

**Our tests have two differentiated lines.**

Control line - green color  
Positive line - red color

Possibility to add **positive and negative controls.**









## Fecal Antigens

### GASTROINTESTINAL VIRUS

#### ADENOVIRUS

Reference	Product name	
AD820001V	<b>Adenovirus</b> (20 test cassettes with stool vials)	
AD820001VC	<b>Adenovirus</b> (20 test cassettes with stool vials + Positive Control)	
AD220001P	<b>Adenovirus</b> 20 Positive Controls Kit	



#### ASTROVIRUS

Reference	Product name	
AT820001V	<b>Astrovirus</b> (20 test cassettes with stool vials)	
AT820001VC	<b>Astrovirus</b> (20 test cassettes with stool vials + Positive Control)	
AT220001P	<b>Astrovirus</b> 20 Positive Controls Kit	



#### ENTEROVIRUS

Reference	Product name	
EV820001V	<b>Enterovirus</b> (20 test cassettes with stool vials)	
EV820001VC	<b>Enterovirus</b> (20 test cassettes with stool vials + Positive Control)	
EV220001P	<b>Enterovirus</b> 20 Positive Controls Kit	







#### NOROVIRUS

Reference	Product name	
NN820001V	<b>Norovirus</b> (20 test cassettes with stool vials)	
NN820001VC	<b>Norovirus</b> (20 test cassettes with stool vials + Positive Control)	
NP220001P	<b>Norovirus</b> 20 Positive Controls Kit	





 **ROTAVIRUS**

Reference	Product name	
R820001V	<b>Rotavirus</b> (20 test cassettes with stool vials)	
R820001VC	<b>Rotavirus</b> (20 test cassettes with stool vials + Positive Control)	
R220001P	<b>Rotavirus</b> 20 Positive Controls Kit	



 **COMBO**

Reference	Product name	
XR882001V	<b>Rotavirus + Adenovirus</b> (20 test combo cassettes with stool vials)	
XR882001VC	<b>Rotavirus + Adenovirus</b> (20 test combo cassettes with stool vials + Positive Controls)	
NV882001V	<b>Norovirus GI + GII</b> (20 test combo cassettes with stool vials)	
NV882001VC	<b>Norovirus GI + GII</b> (20 test combo cassettes with stool vials + Positive Controls)	
RN882001V	<b>Rotavirus + Noro GI/GII</b> (20 test combo cassettes with stool vials)	
RN882001VC	<b>Rotavirus + Noro GI/GII</b> (20 test combo cassettes with stool vials + Positive Controls)	

 **TRIPLE**

Reference	Product name	
RT872001V	<b>Rotavirus + Adenovirus + Astrovirus</b> (20 test triple cassettes with stool vials)	
RT872001VC	<b>Rotavirus + Adenovirus + Astrovirus</b> (20 test triple cassettes with stool vials + Positive Controls)	
RA872001V	<b>Rotavirus + Adenovirus + Norovirus</b> (20 test triple cassettes with stool vials)	
RA872001VC	<b>Rotavirus + Adenovirus + Norovirus</b> (20 test triple cassettes with stool vials + Positive Controls)	

 **QUAD**



Reference	Product name	
RG862001V	<b>Rota + Adeno + Astro + Noro</b> (20 test quad cassettes with stool vials)	
RG862001VC	<b>Rota + Adeno + Astro + Noro</b> (20 test quad cassettes with stool vials + Positive Controls)	





## Fecal Antigens

### GASTROINTESTINAL BACTERIA

#### HELICOBACTER PYLORI

Reference	Product name	
P820001V	<b>H. Pylori</b> (20 test cassettes with stool vials)	
P820001VC	<b>H. Pylori</b> (20 test cassettes with stool vials + Positive Control)	
P220001P	<b>H. Pylori</b> 20 Positive Controls Kit	



#### CLOSTRIDIUM DIFFICILE

Reference	Product name	
GD820001V	<b>Clostridium Difficile antigen GDH</b> (20 test cassettes with stool vials)	
GD820001VC	<b>Clostridium Difficile antigen GDH</b> (20 test cassettes with stool vials + Positive Control)	
GD220001P	<b>Clostridium Difficile antigen GDH</b> 20 Positive Controls Kit	
TA220001P	<b>Clostridium Difficile Toxin A</b> 20 Positive Controls Kit	
TB220001P	<b>Clostridium Difficile Toxin B</b> 20 Positive Controls Kit	

#### ESCHERICHIA COLI

Reference	Product name	
E820001V	<b>E. Coli O157:H7</b> (20 test cassettes with stool vials)	
E820001VC	<b>E. Coli O157:H7</b> (20 test cassettes with stool vials + Positive Control)	
E220001P	<b>E. Coli O157:H7</b> 20 Positive Controls Kit	
VT220001P	<b>EHEC VT1</b> Verotoxin 1, 20 Positive Controls Kit	
VX220001P	<b>EHEC VT2</b> Verotoxin 2, 20 Positive Controls Kit	

#### CAMPYLOBACTER

Reference	Product name	
CA820001V	<b>Campylobacter</b> (20 test cassettes with stool vials)	
CA820001VC	<b>Campylobacter</b> (20 test cassettes with stool vials + Positive Control)	
CA220001P	<b>Campylobacter</b> 20 Positive Controls Kit	

## SALMONELLA

Reference	Product name
SA820001V	<b>Salmonella</b> (20 test cassettes with stool vials)
SA820001VC	<b>Salmonella</b> (20 test cassettes with stool vials + Positive Control)
SA220001P	<b>Salmonella</b> 20 Positive Controls Kit

## SALMONELLA TYPHI

Reference	Product name
ST820001V	<b>Salmonella Typhi</b> (20 test cassettes with stool vials)
ST820001VC	<b>Salmonella Typhi</b> (20 test cassettes with stool vials + Positive Control)
ST220001P	<b>Salmonella Typhi</b> 20 Positive Controls Kit

## SALMONELLA PARATYPHI

Reference	Product name
SP820001V	<b>Salmonella Paratyphi</b> (20 test cassettes with stool vials)
SP820001VC	<b>Salmonella Paratyphi</b> (20 test cassettes with stool vials + Positive Control)
SP220001P	<b>Salmonella Paratyphi</b> 20 Positive Controls Kit

## YERSINISA O:3

Reference	Product name
YE820001V	<b>Yersinia Enterocolitica O:3</b> (20 test cassettes with stool vials)
YE820001VC	<b>Yersinia Enterocolitica O:3</b> (20 test cassettes with stool vials + Positive Control)
YE220001P	<b>Yersinia Enterocolitica O:3</b> 20 Positive Controls Kit

## YERSINISA O:9



Reference	Product name
YC820001V	<b>Yersinia Enterocolitica O:9</b> (20 test cassettes with stool vials)
YC820001VC	<b>Yersinia Enterocolitica O:9</b> (20 test cassettes with stool vials + Positive Control)
YC220001P	<b>Yersinia Enterocolitica O:9</b> 20 Positive Controls Kit



## Fecal Antigens




### COMBO

Reference	Product name	
CD882001V	<b>Clostridium Difficile Toxin A + Toxin B</b> (20 test combo cassettes with stool vials)	
CD882001VC	<b>Clostridium Difficile Toxin A + Toxin B</b> (20 test combo cassettes with stool vials + Positive Controls)	
PT882001VF	<b>H. Pylori + Transferrin</b> (20 test combo cassettes with stool vials)	
PT882001VFC	<b>H. Pylori + Transferrin</b> (20 test combo cassettes with stool vials + Positive Controls)	
GF882001VF	<b>GDH + Lactoferrin</b> (20 test combo cassettes with stool vials)	
GF882001VFC	<b>GDH + Lactoferrin</b> (20 test combo cassettes with stool vials + Positive Controls)	
SX882001V	<b>EHEC VT1 + VT2</b> (20 test combo cassettes with stool vials)	
SX882001VC	<b>EHEC VT1 + VT2</b> (20 test combo cassettes with stool vials + Positive Controls)	
TP882001V	<b>Salmonella Typhi + Paratyphi</b> (20 test combo cassettes with stool vials)	
TP882001VC	<b>Salmonella Typhi + Paratyphi</b> (20 test combo cassettes with stool vials + Positive Controls)	
YS882001V	<b>Yersinia Enterocolitica O:3 + O:9</b> (20 test combo cassettes with stool vials)	
YS882001VC	<b>Yersinia Enterocolitica O:3 + O:9</b> (20 test combo cassettes with stool vials + Positive Controls)	



### TRIPLE

Reference	Product name	
GX872001V	<b>C. Difficile GDH + Toxin A + Toxin B</b> (20 test triple cassettes with stool vials)	
GX872001VC	<b>C. Difficile GDH + Toxin A + Toxin B</b> (20 test triple cassettes with stool vials + Positive Controls)	
EX872001V	<b>E. Coli O157:H7 + VT1 + VT2</b> (20 test triple cassettes with stool vials)	
EX872001VC	<b>E. Coli O157:H7 + VT1 + VT2</b> (20 test triple cassettes with stool vials + Positive Controls)	





### QUAD



Reference	Product name
LX862001VF	<b>C. Difficile GDH + Toxin A + Toxin B + Lactoferrin</b> (20 test quad cassettes with stool vials)
LX862001VFC	<b>C. Difficile GDH + Toxin A + Toxin B + Lactoferrin</b> (20 test quad cassettes with stool vials + Positive Controls)

## GASTROINTESTINAL PARASITES



### CRYPROSPORIDIUM

Reference	Product name	
K820001V	<b>Crypto</b> (20 test cassettes with stool vials)	
K820001VC	<b>Crypto</b> (20 test cassettes with stool vials + Positive Control)	
K220001P	<b>Crypto</b> 20 Positive Controls Kit	



### GIARDIA

Reference	Product name	
GW820001V	<b>Giardia+</b> (20 test cassettes with stool vials)	
GW820001VC	<b>Giardia+</b> (20 test cassettes with stool vials + Positive Control)	
GW220001P	<b>Giardia+</b> 20 Positive Controls Kit	



### ENTAMOEBIA

Reference	Product name	
EH820001V	<b>Entamoeba</b> (20 test cassettes with stool vials)	
EH820001VC	<b>Entamoeba</b> (20 test cassettes with stool vials + Positive Control)	
EH220001P	<b>Entamoeba</b> 20 Positive Controls Kit	

### COMBO

Reference	Product name	
NW882001V	<b>Crypto + Giardia</b> (20 test combo cassettes with stool vials)	
NW882001VC	<b>Crypto + Giardia</b> (20 test combo cassettes with stool vials + Positive Controls)	

### TRIPLE

Reference	Product name	
NE872001V	<b>Crypto + Giardia + Entamoeba</b> (20 test combo cassettes with stool vials)	
NE872001VC	<b>Crypto + Giardia + Entamoeba</b> (20 test combo cassettes with stool vials + Positive Controls)	



## Respiratory Antigens

### RESPIRATORY VIRUS

#### INFLUENZA

Reference	Product name
Y820001PC	<b>Influenza A</b> (20 test cassettes with stool vials + Positive Control)
Y220001P	<b>Influenza A</b> 20 Positive Controls Kit
YB220001P	<b>Influenza B</b> 20 Positive Controls Kit

#### RSV (RESPIRATORY SYNCYTIAL VIRUS)

Reference	Product name
RV820001PC	<b>RSV</b> (20 test cassettes with stool vials + Positive Control)
RV220001P	<b>RSV</b> 20 Positive Controls Kit

#### ADENOVIRUS RESPIRATORY

Reference	Product name
AP820001PC	<b>Adenovirus Resp.</b> (20 test cassettes with stool vials + Positive Control)
AP220001P	<b>Adenovirus Resp.</b> 20 Positive Controls Kit

#### SARS COV-2

Reference	Product name
SC820001P	<b>SARS CoV-2</b> (20 test cassettes with stool vials)
SC820001PC	<b>SARS CoV-2</b> (20 test cassettes with stool vials + Positive Control)
SC860001PC	<b>SARS CoV-2</b> (60 test cassettes with stool vials + Positive Control)
SC220001P	<b>SARS CoV-2</b> 20 Positive Controls Kit



## COMBO

Reference	Product name
Z882001PC	<b>Influenza A + B</b> (20 test cassettes with stool vials + Positive Control)
VR882001PC	<b>RSV + Adenovirus Resp.</b> (20 test cassettes with stool vials + Positive Control)



## TRIPLE

Reference	Product name
ZS872001PC	<b>Influenza A + B + RSV</b> (20 test cassettes with stool vials + Positive Control)
SS872001PC	<b>SARS CoV-2 + Flu A + Flu B</b> (20 test cassettes with stool vials + Positive Control)



## QUAD

Reference	Product name
ZV862001PC	<b>Influenza A + B + RSV + Adeno Resp.</b> (20 test cassettes with stool vials + Positive Control)
SV862001PC	<b>SARS CoV-2 + Flu A + Flu B + RSV</b> (20 test cassettes with stool vials + Positive Control)

## RESPIRATORY BACTERIA

### STREP. A

Reference	Product name
T820201PC	<b>Strep A</b> (20 test cassettes with diluent bottles + tubes and pipettes vials + Positive Control)
T220001P	<b>Strep A</b> 20 Positive Controls Kit

Any test KIT which includes positive control could have negative control too adding C at the end of the reference.

Reference	Product name
CN320001P	Negative Control 20 Negative Controls Kit



## Urine Antigens

### URINE BACTERIA

#### LEGIONELLA PNEUMOPHILA

Reference	Product name
LP820201PC	<b>Legionella</b> (20 test cassettes with diluent bottle + Positive Control + Negative Control)
LP220001P	<b>Legionella</b> 20 Positive Controls Kit

#### STREPTOCOCCUS PNEUMONIAE

Reference	Product name
SN820201PC	<b>S. pneumoniae</b> (20 test cassettes with diluent bottle + Positive Control + Negative Control)
SN220001P	<b>S. pneumoniae</b> 20 Positive Controls Kit

#### COMBO

Reference	Product name
NL882201PC	<b>S. pneumoniae + Legionella</b> (20 test cassettes with diluent bottle + Positive Controls + N. C.)



## Inflammatory & Tumor Makers

#### FOB

Reference	Product name
F820001F	<b>FOB</b> (20 test cassettes with stool vials)
F820001FC	<b>FOB</b> (20 test cassettes with stool vials + Positive Control)
F220001P	<b>FOB</b> 20 Positive Controls Kit

#### CALPROTECTIN

Reference	Product name
CP820001F	<b>Calprotectin</b> (20 test cassettes with stool vials)
CP820001FC	<b>Calprotectin</b> (20 test cassettes with stool vials + Positive Control)
CP220001P	<b>Calprotectin</b> 20 Positive Controls Kit

 **TRANSFERRIN**

Reference	Product name
TF820001F	<b>Transferrin</b> (20 test cassettes with stool vials)
TF820001FC	<b>Transferrin</b> (20 test cassettes with stool vials + Positive Control)
TF220001P	<b>Transferrin</b> 20 Positive Controls Kit

 **LACTOFERRIN**

Reference	Product name
LT820001F	<b>Lactoferrin</b> (20 test cassettes with stool vials)
LT820001FC	<b>Lactoferrin</b> (20 test cassettes with stool vials + Positive Control)
LT220001P	<b>Lactoferrin</b> 20 Positive Controls Kit

 **COMBO**

Reference	Product name
FB882001F	<b>FOB 50 + 200</b> (20 test combo cassettes with stool vials)
FB882001FC	<b>FOB 50 + 200</b> (20 test combo cassettes with stool vials + Positive Controls)
CC882001F	<b>Calprotectin 50 + 200</b> (20 test combo cassettes with stool vials)
CC882001FC	<b>Calprotectin 50 + 200</b> (20 test combo cassettes with stool vials + Positive Controls)
FT882001F	<b>FOB + Transferrin</b> (20 test combo cassettes with stool vials)
FT882001FC	<b>FOB + Transferrin</b> (20 test combo cassettes with stool vials + Positive Controls)
CL882001F	<b>Calprotectin + Lactoferin</b> (20 test combo cassettes with stool vials)
CL882001FC	<b>Calprotectin + Lactoferin</b> (20 test combo cassettes with stool vials + Positive Controls)

 **QUAD**

Reference	Product name
FC862001F	<b>FOB + Transferrin + Calprotectin + Lactoferrin</b> (20 test combo cassettes with stool vials)
FC862001FC	<b>FOB + Transferrin + Calprotectin + Lactoferrin</b> (20 test combo cassettes with stool vials + Positive Controls)

**Certest Immunodiagnostic** offers simple, sensitive, specific screening assays for a presumptive diagnosis of infectious diseases.

## Advantages:



Noninvasive diagnostic



All included. No additional equipment needed



Low operational costs



Immediate results



Ease of use and interpretation

From Discovery  
to Delivery,  
**Biotech** you  
can trust.

**certest**

IMMUNODIAGNOSTIC

**Certest Biotec, S.L.**

Pol. Industrial Río Gállego II · Calle J, N°1  
50840, San Mateo de Gállego, Zaragoza (Spain)  
Tel. (+34) 976 520 354  
salesidx@certest.es  
[www.certest.es](http://www.certest.es)



Find us on

@certestbiotec

Certest Rapid Test/GEN-0426EN

Modification rights reserved. All rights reserved. © Certest Biotec, S.L.  
The products, services and data set out in this document may suffer changes and/or variations on the texts and pictures shown.