



TECHNICAL SHEET

PRODUCT NAME: SARS-CoV-2 recombinant Spike Glycoprotein (S) (full sequence) (mammalian expression)

Code: MT-25C19S

Physical State: Liquid

Source: *Hek293 cells*

Description: SARS-CoV-2 recombinant Spike Glycoprotein (S) (full sequence) (mammalian expression) expressed and purified from *Hek293 cells*.

Appearance: Clear

Purity: >90%

Storage Conditions: Short term storage at 2-8°C. Long term storage at -20°C. Avoid multiple/freeze thaw cycles by storing multiple aliquots at -20°C.

Shipping: With cool packs.

Health & Safety Information: Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.

Presentation: Tris buffer with additives and salts, pH: 8.00, 0.09% Sodium Azide (NaN₃).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 19/08/2022