

Helicobacter stool strips

- CLIA technique is more sensitive than other immunological techniques, which is a great advantage in the rapid and accurate diagnosis of infectious diseases.

Helicobacter pylori is a gram-negative bacteria that cause stomach infections. It is found in nearly two thirds of the world populations. It is involved in peptic ulcers, and in 80% of the stomach cancer. Infection usually happened in the childhood, but pathology is revealed in adulthood.

Stool antigen detection of *H. pylori* is becoming a high demanded non-invasive diagnostic test. CerTest CLIA Pylori are the best choice to complete your diagnosis panel.

CerTest CLIA strip offers:



- Broad panel for stool antigen detection.
- Ready to use format, minimum handling needed.
- Adaptable in open CLIA plate readers.





Helicobacter pylori stool strips

►► Specifications:

Technique	Chemiluminescence Immunoassay (HRP)
Assay format	Enzyme chemiluminescence immunoassay
Test Type	Qualitative
Samples	Faecal
Number of determinations	60 (minimum)
Storage	2-8°C
Reference method	ELISA ProSpect Campylobacter micrplate (OXOID)
Clinical sensitivity	100%
Clinical specificity	97,8%

►► Ordering Info:

Ref	Name	Content
PYL-Ag	H. pylori Ag Plate	<ul style="list-style-type: none">• Coated strips (12x)• HRP-conjugated antibody (1x)• Positive control (2x)• Negative control (2x)• Cut off (2x)• Substrate A (1x)• Substrate B (1x)
WASH-3-100mL	Wash buffer	<ul style="list-style-type: none">• 100 mL 20X concentrated buffer
VFE-002	CerTest CLIA H. pylori Faecal Sample Vials kit	<ul style="list-style-type: none">• Faecal extraction vials (100 vials)

►► Compatibility with:

ThunderBolt® and Bolt™ (Gold Standard Diagnostics)

ChemWell Fusion (Awareness Technology)

Available in BULK format 

For more information and use procedure, read the instructions for use included in this product.

CerTest
BIOTEC

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
Tel. (+34) 976 520 354 | clia@certest.es

www.certest.es



CLIA/PYL/Plate/0622EN