

CerTest CLIA H. pylori Ag

►► Integral solution

H. pylori Ag is a stool antigen detection kit. Sensitivity and specificity for early diagnosis with CerTest CLIA.

CLIA technique is more sensitive than other immunological techniques, which is a great advantage in the rapid and accurate diagnosis of infectious diseases.

Helicobacter pylori is a gram-negative bacteria that cause stomach infections. It is found in nearly two thirds of the world populations. It is involved in peptic ulcers, and in 80% of the stomach cancer. Infection usually happened in the childhood, but pathology is revealed in adulthood.

Stool antigen detection of *H. pylori* is becoming a high demanded non-invasive diagnostic test. CerTest CLIA Pylori are the best choice to complete your diagnosis panel.



CerTest CLIA Integran solution offers:



- **Stool, respiratory and serological panels** for infectious diseases diagnostics
- **Ready to use reagents**, minimum handling needed
- **Batch and random a Sample access** available





CerTest CLIA H. pylori Ag

►► Specifications:

Technique	Chemiluminescence Immunoassay (AE)
Assay format	Single sandwich format
Test Type	Qualitative
Samples	Faecal
Platform	Shine i1000
Number of determinations	200
Storage	2-8°C
Reference method	Amplified IDEIA Hp Star (Oxoid)
Clinical sensitivity	100%
Clinical specificity	96%

►► Ordering Info:

Ref	Name	Content
PYL001	CerTest CLIA H. pylori Ag 001	<ul style="list-style-type: none">H. pylori Ag Beads & Conjugate S-cartridge2 level controls
PYLCal	CerTest CLIA H. pylori Ag Calibration set	<ul style="list-style-type: none">1 level calibrator
CL007	CLIA Analyzer	<ul style="list-style-type: none">Model Shine i1910
VFE-002	CerTest CLIA H. pylori Faecal Sample Vials kit	<ul style="list-style-type: none">Fecal extraction vials (100 vials)

Available in BULK format 

For more information and use procedure, read the instructions for use included in this product.

CerTest
BIOTEC

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
Tel. (+34) 976 520 354 | clia@certest.es

www.certest.es



CLIA/PYL/IS/1122EN