


PLUS C-Reactive Protein Turbilatex® Combo

**C-reactive protein (CRP)
is a suitable marker
of inflammation process**

thanks to its long half-life,
independence from age,
long stability and the availability
of a wide range of number
of diagnostic methods available.



CerTest
BIOTEC

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, N°1
50840, San Mateo de Gállego, Zaragoza (Spain)

Tel. (+34) 976 520 354

Fax (+34) 976 106 268

certest@certest.es | turbilatex@certest.es

www.certest.es

PLUS C-Reactive Protein Turbilatex® Combo

This assay is indicated to evaluate the amount of C-reactive protein in serum samples.

Detects HS-CRP and CRP in a single kit, contains the components to realize both parameters.

►► HIGH SENSITIVITY C-REACTIVE PROTEIN

Cardiovascular accidents are the main cause of death in developed countries. It has been established that inflammation is a fundamental process in the pathogenesis of arteriosclerosis in all its stages, initiation, growth, and complication of arteriosclerotic plaque, considered an inflammatory response to the damage.


High-sensitivity C-reactive protein (HS-CRP) concentration increases the risk of myocardial infarction and cerebrovascular disease.

►► C-REACTIVE PROTEIN

C-reactive protein (CRP) participates in innate immunity as a first defense system against infection.

It is an early indicator of infection or inflammation because its levels begin to rise in the blood shortly after the onset of infection or inflammation before the onset of other symptoms such as fever or pain and can increase up to 50,000-fold in acute inflammatory states.

PLUS CRP
human serum samples



References	Description
TL-022CRP200ED	Plus C-Reactive Protein Turbilatex® Combo kit 200 det.
TL-022CRP800ED	Plus C-Reactive Protein Turbilatex® Combo kit 800 det.
TL-022HS200ED	High Sensitivity C-Reactive Protein Turbilatex® Combo kit 200 det.
TL-022HS800ED	High Sensitivity C-Reactive Protein Turbilatex® Combo kit 800 det.
TL-022CR200ED	C-Reactive Protein Turbilatex® Combo kit 200 det.
TL-022CR800ED	C-Reactive Protein Turbilatex® Combo kit 800 det.

- PLUS CRP Turbilatex®, High Sensitivity CRP Turbilatex®, CRP Turbilatex® are a latex turbidimetric assays only for the quantitative detection of C-reactive protein in human serum samples (not to be used for body fluid as whole blood or plasma).
- These assays are simple and widely applicable. These products are optimized for several automated analyser.
- For professional in vitro diagnostic use only.