



TECHNICAL SHEET

PRODUCT NAME: Anti-IMP mAb (clone IM01)

Code:	MT-16IM01, MT-16IM01s1
Clone:	IM01
Isotype:	IgG ₁
Physical State:	Liquid
Host Animal (source):	Mouse
Specificity:	IM01 is specific for IMP carbapenemase detection.
Cross Reactivity:	No cross reactivity with other carbapenemases.
Purification:	Anti-IMP mAb (clone IM01) has been purified in a multi-step procedure. ≥95% pure by SDS-PAGE.
Storage Conditions:	Storage at 2-8°C.
Health & Safety Information:	This product contains sodium azide as preservative. Although the amount is 0.09%, this substance is poisonous and hazardous and it should be handled by trained staff only. Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.
Presentation:	PBS pH: 7.4, 0.09% Sodium Azide (NaN ₃).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 01/06/2022