



TECHNICAL SHEET

PRODUCT NAME: Anti-VIM mAb (clone VI21)

Code:	MT-16VI21, MT-16VI21s1
Clone:	VI21
Isotype:	IgG ₁
Physical State:	Liquid
Host Animal (source):	Mouse
Specificity:	VI21 is specific for VIM carbapenemase detection.
Cross Reactivity:	No cross reactivity with other carbapenemase.
Purification:	Anti-VIM mAb (clone VI21) has been purified in a multi-step procedure. ≥95% pure by SDS-PAGE.
Storage Conditions:	Storage at 2-8°C.
Health & Safety Information:	This product contains sodium azide as preservative. Although the amount is 0.09%, this substance is poisonous and hazardous and it should be handled by trained staff only. Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.
Presentation:	PBS pH: 7.4, 0.09% Sodium Azide (NaN ₃).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 01/06/2022