

VIASURE

RT-PCR Complete Solution

V-Lab96

Open platform for *in vitro* diagnostics

VIASURE V-Lab96 allows users to analyze 96 samples simultaneously for qualitative and quantitative Real Time PCR.



Viasure
Real Time
PCR platform



Improved optical system for **high specificity and sensitive** detection.



Reduced operating time.



High sample throughput in diagnostic lab applications.



Data registering and accounting in electronic format for LIS connexion.



Power-off protection. Avoid the loss of experimental data and waste of reagents due to abrupt power-off.



Analysis performed in a closed test tube to minimize possible contamination of further experiments.



Integrated 10.4" **TouchScreen.**



Compatible with all VIASURE Real Time PCR Detection Kits.



Qualitative & Quantitative data analysis.



Automatic interpretation and analysis of results



VIASURE V-Lab96
software interface

