



TECHNICAL SHEET

PRODUCT NAME: Anti-NDM mAb (clone ND24)

Code:	MT-16ND24, MT-16ND24s1
Clone:	ND24
Isotype:	IgG ₁
Physical State:	Liquid
Host Animal (source):	Mouse
Specificity:	ND24 is specific for NDM carbapenemase detection.
Cross Reactivity:	No cross reactivity with others carbapenemases.
Purification:	Anti-NDM mAb (clone ND24) has been purified in a multi-step procedure. ≥95% pure by SDS-PAGE.
Storage Conditions:	Storage at 2-8°C.
Health & Safety Information:	This product contains sodium azide as preservative. Although the amount is 0.09%, this substance is poisonous and hazardous and it should be handled by trained staff only. Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.
Presentation:	PBS pH: 7.4, 0.09% Sodium Azide (NaN ₃).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 04/07/2022