

Technical Sheet

Product Name: Anti-IMP mAb (clone IMP18)

Business Unit: CerTest BioSCIENCE

TECHNICAL DETAILS:

Code: MT-16IM18

Physical State: Liquid

Source: Mice

Description: Anti-IMP mAb (clone IMP18)

Appearance: Clear

Purity: > 95% (determined by SDS gel electrophoresis)

Storage Conditions: Shipping: 2-8°C. Short term storage at 2-8°C. Long term storage at -20°C. Avoid multiple/freeze thaw cycles by storing multiple aliquots at -20°C.

Health & Safety Information: This product contains sodium azide as preservative. Although the amount is 0.09%, this substance is poisonous and hazardous and it should be handled by trained staff only. Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.

Presentation: PBS 1X, pH 7.4, 0.09% Sodium Azide (NaN₃).

Purification Method: Affinity chromatography (Protein A)


Isotype: IgG1

Molecular Weight (MW): 150KDa


Isoelectric point (pI): 6.5-7

Denaturation Midpoint Temperature (TM): 73.7±0.6 °C

SDS IMAGE



HPLC IMAGE



Date of Technical sheet: **07/12/2022**

Signed:



Rafael Clavería
Bioscience Manager

Code of the document and version: **Tshb-MT-16IM18 v.00**

CHANGE CONTROL:

| CHANGE CONTROL | | |
|----------------|------------------------------|------------|
| Version nº | Changes | Date |
| 00 | Original version of the file | 07/12/2022 |