

Technical Sheet

Product Name:	MCR-1 recombinant protein
Business Unit:	CerTest Raw Material

TECHNICAL DETAILS:

Code:	MT-25MCR1, MT-25MCR1s1
Physical State:	Liquid
Source:	<i>Escherichia coli</i>
Description:	MCR-1 recombinant protein expressed and purified from <i>Escherichia coli</i> .
Storage Conditions:	Shipping: 2-8°C. Short term storage at 2-8°C. Long term storage at -20°C. Avoid multiple/freeze thaw cycles by storing multiple aliquots at -20°C
Purity:	>90%
Health & Safety Information:	This product contains sodium azide as preservative. Although the amount is 0.09%, this substance is poisonous and hazardous and it should be handled by trained staff only. Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations
Presentation:	Tris buffer with additives and salts, pH: 8.0, 0.09% Sodium Azide (NaN ₃).

Date of Technical sheet: **04/08/2023**

Signed:



Rafael Clavería
Raw Material Manager

Code of the document and version: **Tshb-MT-25MCR1, MT-25MCR1s1 v.00**

CHANGE CONTROL:

CHANGE CONTROL		
Version nº	Changes	Date
00	Original version of the file	04/08/2023