


Antibodies & Antigens for immunodiagnostic






certest

RAW MATERIALS

Our expertise in **protein expression & purification** and **monoclonal antibodies development & manufacturing** offer to our customers a priceless support in their challenging projects for IVD product development.

All our products have been properly characterized in our internal platforms and our IVD products, helping us to provide our customer some useful advices about their use and performance.

Additionally, we also offer to **third parties manufacturing** and **R&D services** for the production or the development of your own materials carried out by our high qualified staff in a competitive time.



Index

Antibodies & Antigens for immunodiagnostic

1. Infectious diseases

1.1 Infectious diseases / Respiratory

- 5 1. Strep A
- 5 2. Respiratory Syncytial Virus (RSV)
- 5 3. Adenovirus
- 6 4. SARS-CoV-2 & other Coronavirus
- 7 5. Influenza
- 8 6. Legionella pneumophila
- 8 7. Streptococcus pneumoniae
- 8 8. Mycoplasma pneumoniae
- 8 9. Human metapneumovirus

1.2 Infectious diseases / Gastrointestinal

- 9 1. Campylobacter
- 9 2. Norovirus
- 10 3. H. pylori
- 10 4. Astrovirus
- 10 5. Giardia intestinalis
- 10 6. Cryptosporidium parvum
- 11 7. Clostridium difficile
- 11 8. Rotavirus
- 12 9. Salmonella
- 12 10. Listeria
- 12 11. Enterovirus
- 13 12. Entamoeba histolytica/dispar
- 13 13. Yersinia enterocolitica
- 13 14. Shigella
- 13 15. E. coli O157

2. Antimicrobial resistance


- 14 1. MCR-1
- 14 2. OXA-48
- 14 3. KPC
- 14 4. NDM
- 14 5. VIM
- 15 6. IMP

3. Tumor & Inflammation markers

- 15 1. Human Calprotectin
- 15 2. Human Haemoglobin (FOB)
- 16 3. Human Pancreatic Elastase
- 16 4. Alfa-1-antitrypsin recombinant protein
- 16 5. Human Transferrin
- 16 6. Human IL-6
- 17 7. Human Lactoferrin
- 17 8. Human S100b
- 17 9. Human S100A12
- 17 10. Human S100A9 and S100A8

4. Other markers


- 17 1. Ferritin



**Customized developments and productions
for all our clients.**



Features:



Icons:

- | | | | |
|--|---|--|--|
|  Monoclonal Antibodies |  Native Antigens |  Self Pairing | |
|  Polyclonal Antibodies |  VLPs | |  Paired with... |
|  Recombinant Proteins |  Virus | | |
|  Native Antigens Purified | | | |

1.1 Infectious diseases / Respiratory

1. For Strep A detection

 	MT-20TSS	Anti-Strep A pAb (x1mg) / Self Pairing!
	MT-28SAGU	Inactivated Strep A antigen (native extract) (x1mL)

2. For Respiratory Syncytial Virus (RSV) detection

  	MT-16RV11	Anti-RSV mAb (clone RV11) (x1mg) / Self Pairing! /Paired with MT-16RV12 & MT-16RV46
 	MT-16RV12	Anti-RSV mAb (clone RV12) (x1mg) / Paired with MT-16RV11 & MT-16RV49
 	MT-16RV46	Anti-RSV mAb (clone RV46) (x1mg) / Paired with MT-16RV11 & MT-16RV49
 	MT-16RV49	Anti-RSV mAb (clone RV49) (x1mg) / Paired with MT-16RV12 & MT-16RV46
	MT-25RSV	RSV recombinant fusion protein (x1mg)
	MT-29RVV	Inactivated RSV antigen (native extract) (x1mL)

3. For Adenovirus detection

 	MT-16A68	Anti-Adenovirus mAb (clone A68) (x1mg) / Paired with MT-16A93
 	MT-16A93	Anti-Adenovirus mAb (clone A93) (x1mg) / Paired with MT-16A68 & MT-16A114
  	MT-16A114	Anti-Adenovirus mAb (clone A114) (x1mg) / Self Pairing! / Paired with MT-16A93 & MT-16A115
 	MT-16A115	Anti-Adenovirus mAb (clone A115) (x1mg) / Paired with MT-16A114
	MT-25HEX	Adenovirus HEXON protein (x1mg)
	MT-29ADU	Inactivated Adenovirus antigen (native extract) (x1mL)

 Monoclonal Antibodies

 Polyclonal Antibodies

 Recombinant Proteins

 Native Antigens

 Virus

 Paired with...

 Self Pairing

1.1 Infectious diseases / Respiratory

4. For SARS-CoV-2 detection

4.1 Monoclonal antibodies for SARS-CoV-2 detection (Nucleoprotein detection)

		MT-16CV01	Anti SARS-CoV-2 mAb (clone CV01) (x1mg) / Paired with MT-16CV40 & MT-16CV74
		MT-16CV10	Anti SARS-CoV-2 mAb (clone CV10) (x1mg) / Paired with MT-16CV40 & MT-16CV74
		MT-16CV40	Anti SARS-CoV-2 mAb (clone CV40) (x1mg) / Paired with MT-16CV01 & MT-16CV10
		MT-16CV74	Anti SARS-CoV-2 mAb (clone CV74) (x1mg) / Paired with MT-16CV01 & MT-16CV10

4.2 Monoclonal antibodies for SARS-CoV-2 detection (Spike protein & RBD detection)

		MT-16SV53	Anti SARS-CoV-2 mAb (clone SV53) (x1mg) / Self Pairing!
---	---	-----------	---

4.3 Antigens. Mammalian expression

	MT-25C19NPm	SARS-CoV-2 recombinant Nucleoprotein (N) (full sequence) (x1mg)
	MT-25C19S	SARS-CoV-2 recombinant Spike Glycoprotein (S) (full sequence) (x1mg)
	MT-25C19S.b.1.1.7	SARS-CoV-2 recombinant Spike Glycoprotein (S) Alpha variant (x1mg)
	MT-25C19S.b.1.351	SARS-CoV-2 recombinant Spike Glycoprotein (S) Beta variant (x1mg)
	MT-25C19S.p.1	SARS-CoV-2 recombinant Spike Glycoprotein (S) Gamma variant (x1mg)
	MT-25C19S.b.1.617.2	SARS-CoV-2 recombinant Spike Glycoprotein (S) Delta variant (x1mg)
	MT-25C19S.BA.1	SARS-CoV-2 recombinant Spike Glycoprotein (S) Omicron BA.1 variant (x1mg)
	MT-25C19S.BA.2	SARS-CoV-2 recombinant Spike Glycoprotein (S) Omicron BA.2 variant (x1mg)
	MT-25RBD	SARS-CoV-2 recombinant Receptor Binding Domain (RBD) (mammalian expression) (x1mg)



All known SARS-CoV-2 variants are
recognized by our antibodies!



Monoclonal Antibodies



Recombinant Proteins



Paired with...



Self Pairing

4.4 Bacterial expression

	MT-25C19NP	SARS-CoV-2 recombinant Nucleoprotein (N) (full sequence) (x1mg)
	MT-25C19NP.b.1.1.7	SARS-CoV-2 recombinant Nucleoprotein (N) Alpha variant (x1mg)
	MT-25C19NP.b.1.351	SARS-CoV-2 recombinant Nucleoprotein (N) Beta variant (x1mg)
	MT-25C19NP.p.1	SARS-CoV-2 recombinant Nucleoprotein (N) Gamma variant (x1mg)
	MT-25C19NP.b.1.617.2	SARS-CoV-2 recombinant Nucleoprotein (N) Delta variant (x1mg)
	MT-25C19NP.BA.1	SARS-CoV-2 recombinant Nucleoprotein (N) Omicron BA.1 variant (x1mg)
	MT-25C19NP.BA.2	SARS-CoV-2 recombinant Nucleoprotein (N) Omicron BA.2 variant (x1mg)

4.5 Antigens from other Coronavirus

	MT-25SANP	SARS Coronavirus recombinant Nucleoprotein (NP) (full sequence) (x1mg)
	MT-25MENP	MERS Coronavirus recombinant Nucleoprotein (NP) (full sequence) (x1mg)
	MT-25229NP	229E Coronavirus recombinant Nucleoprotein (NP) (full sequence) (x1mg)
	MT-25OCNP	OC43 Coronavirus recombinant Nucleoprotein (NP) (full sequence) (x1mg)
	MT-25HKNP	HKU1 Coronavirus recombinant Nucleoprotein (NP) (full sequence) (x1mg)
	MT-25NLNP	NL63 Coronavirus recombinant Nucleoprotein (NP) (full sequence) (x1mg)

5. For Influenza detection

5.1 For Influenza A detection

		MT-18Y77	Anti-Influenza A mAb (clone Y77) (x1mg) / Self Pairing!
		MT-16IA57	Anti-Influenza A mAb (clone IA57) (x1mg) / Paired with MT-16IA61
		MT-16IA61	Anti-Influenza A mAb (clone IA61) (x1mg) / Paired with MT-16IA57
		MT-25FAN	Influenza A recombinant nucleoprotein (x1mg)

5.2 For Influenza B detection

		MT-18YB91	Anti-Influenza B mAb (clone YB91) (x1mg) / Self Pairing!
		MT-16IB17	Anti-Influenza B mAb (clone IB17) (x1mg) / Paired with MT-16IB59 & MT-16IB62
		MT-16IB59	Anti Influenza B mAb (clone IB59) (x1 mg) / Paired with MT-16IB17
		MT-16IB62	Anti Influenza B mAb (clone IB62) (x1 mg) / Paired with MT-16IB17
		MT-25FBN	Influenza B recombinant nucleoprotein (x1mg)

1.1 Infectious diseases / Respiratory

6. For Legionella pneumophila detection

		MT-18LN14	Anti-Legionella pneumophila mAb (clone LN14) (x1mg) / Paired with MT-18LN29
		MT-18LN29	Anti-Legionella pneumophila mAb (clone LN29) (x1mg) / Paired with MT-18LN14
		MT-28LNU	Inactivated Legionella pneumophila antigen (native extract) (x1mL)

7. For Streptococcus pneumoniae detection

		MT-18SN3	Anti-Streptococcus pneumoniae mAb (clone SN3) (x1mg) / Paired with MT-18SN4
		MT-18SN4	Anti-Streptococcus pneumoniae mAb (clone SN4) (x1mg) / Paired with MT-18SN3
		MT-28SPNU	Inactivated Streptococcus pneumoniae antigen (native extract) (x1mL)

8. For Mycoplasma pneumoniae detection

		MT-25MPP	Mycoplasma pneumoniae recombinant protein (x1mg)
--	---	----------	--

9. Human metapneumovirus

		MT-25MPF	Human metapneumovirus (HMPV) recombinant fusion protein (x1mg)
		MT-25MPN	Human metapneumovirus (HMPV) recombinant nucleoprotein (x1mg)
		MT-25MPM	Human metapneumovirus (HMPV) recombinant matrix protein (x1mg)



Monoclonal Antibodies



Recombinant Proteins



Native Antigens



Paired with...



Self Pairing

1.2 Infectious diseases / Gastrointestinal

1. For Campylobacter detection

 	MT-16CA29	Anti-Campylobacter mAb (clone CA29) (x1mg) / Paired with MT-16CA30
  	MT-16CA30	Anti-Campylobacter mAb (clone CA30) (x1mg) / Self Pairing! / Paired with MT-16CA29
	MT-25CEP	Campylobacter jejuni recombinant outer membrane protein (x1mg)
	MT-25CCP	Campylobacter coli recombinant outer membrane protein (x1mg)
	MT-28CJU	Inactivated Campylobacter jejuni antigen (native extract) (x1mL)
	MT-28CCU	Inactivated Campylobacter coli antigen (native extract) (x1mL)

2. For Norovirus detection

2.1 For Norovirus GI detection

 	MT-16NG09	Anti-Norovirus GI mAb (clone NG09) (x1mg) / Paired with MT-16NG39
 	MT-18NG28	Anti-Norovirus GI mAb (clone NG28) (x1mg) / Self Pairing!
 	MT-16NG39	Anti-Norovirus GI mAb (clone NG39) (x1mg) / Paired with MT-16NG09
	MT-25NGI1	Norovirus GI.1 recombinant P domain (x1mg)
	MT-25NGI3	Norovirus GI.3 recombinant P domain (x1mg)
	MT-30NGA	Norovirus GI.1 recombinant VLP (x1mg)

2.2 For Norovirus GII detection

 	MT-16NP23	Anti-Norovirus GII mAb (clone NP23) (x1mg) / Self Pairing!
 	MT-16NP124	Anti-Norovirus GII mAb (clone NP124) (x1mg) / Paired with MT-16NP142 & MT-16NP191
 	MT-16NP142	Anti-Norovirus GII mAb (clone NP142) (x1mg) / Paired with MT-16NP124 & MT-16NP191
 	MT-16NP191	Anti-Norovirus GII mAb (clone NP191) (x1mg) / Paired with MT-16NP124 & MT-16NP142
	MT-25NGII4	Norovirus GII.4 recombinant P domain (x1mg)
	MT-25NGII10	Norovirus GII.10 recombinant P domain (x1mg)
	MT-25NGII17	Norovirus GII.17 recombinant P domain (x1mg)
	MT-30NPA	Norovirus GII.4 recombinant VLP (x1mg)

 _Monoclonal Antibodies  _Recombinant Proteins  _VLPs  Native Antigens

 _Paired with...  _Self Pairing

1.2 Infectious diseases / Gastrointestinal

3. For H. pylori detection

 	MT-16P2	Anti-H. pylori mAb (clone P2) (x1mg) / Self Pairing! / Paired with MT-16P59 & MT-16P62
 	MT-16P59	Anti-H. pylori mAb (clone P59) (x1mg) Paired with MT-16P2 & MT-16P62
 	MT-16P62	Anti-H. pylori mAb (clone P62) (x1mg) Paired with MT-16P2 & MT-16P59
	MT-25PCH	H. pylori recombinant outer membrane protein (x1mg)
	MT-28PECU	Inactivated H. pylori antigen (native extract) (x1mL)

4. For Astrovirus detection

 	MT-16AT54	Anti-Astrovirus mAb (clone AT54) (x1mg) / Paired with MT-16AT74, MT-16AT86 & MT-16AT92
 	MT-16AT74	Anti-Astrovirus mAb (clone AT74) (x1mg) / Paired with MT-16AT54
 	MT-16AT86	Anti-Astrovirus mAb (clone AT86) (x1mg) / Paired with MT-16AT54
 	MT-16AT92	Anti-Astrovirus mAb (clone AT92) (x1mg) / Paired with MT-16AT54
	MT-25AST	Astrovirus capsid recombinant protein (x1mg)

5. Giardia intestinalis

5.1 For Giardia intestinalis trophozoite detection

 	MT-16G18	Anti-Giardia mAb trophozoite protein (clone G18) (x1mg) / Paired with MT-16G22
 	MT-16G22	Anti-Giardia mAb trophozoite protein (clone G22) (x1mg) / Paired with MT-16G18
	MT-25A1G	Giardia intestinalis trophozoite recombinant protein (x1mg)

5.2 For Giardia intestinalis cyst detection

 	MT-16GR7	Anti-Giardia mAb cyst protein (clone GR7) (x1mg) / Paired with MT-16GR16
 	MT-16GR16	Anti-Giardia mAb cyst protein (clone GR16) (x1mg) / Paired with MT-16GR7
	MT-25GCP	Giardia intestinalis cyst recombinant protein (x1mg)

6. For Cryptosporidium parvum detection

 	MT-16CR23	Anti-Crypto mAb (clone CR23) (x1mg) / Self Pairing!
	MT-29KOE	Inactivated Cryptosporidium parvum antigen (native extract) (x1mL)



Monoclonal Antibodies



Recombinant Proteins



Native Antigens



Paired with...



Self Pairing

1.2 Infectious diseases / Gastrointestinal

7. For Clostridium difficile detection

7.1 For Clostridium difficile (GDH) detection

		MT-16GD10	Anti-GDH mAb (clone GD10) (x1mg) / Self Pairing!
		MT-25GDH	Clostridium difficile GDH recombinant protein (x1mg)

7.2 For Clostridium difficile Toxin A detection

		MT-16TA22	Anti-CD Toxin A mAb (clone TA22) (x1mg) / Paired with MT-16TA38	
		MT-16TA35	Anti-CD Toxin A mAb (clone TA35) (x1mg) / Paired with MT-16TA38	
			MT-16TA38	Anti-CD Toxin A mAb (clone TA38) (x1mg) / Self Pairing! / Paired with MT-16TA22 & MT-16TA35
		MT-25TXA	C. difficile Toxin A recombinant protein (x1mg) (fragment without toxic activity)	

7.3 For Clostridium difficile Toxin B detection

		MT-16TB7	Anti-CD Toxin B mAb (clone TB7) (x1mg) / Paired with MT-16TB8
		MT-16TB8	Anti-CD Toxin B mAb (clone TB8) (x1mg) / Paired with MT-16TB7 & MT-16TB75
		MT-18TB41	Anti-CD Toxin B mAb (clone TB41) (x1mg) / Paired with MT-16TB75
		MT-16TB75	Anti-CD Toxin B mAb (clone TB75) (x1mg) / Paired with MT-16TB8, MT-18TB41, MT-16TB175 & MT-16TB194
		MT-16TB175	Anti-CD Toxin B mAb (clone TB175) (x1mg) / Paired with MT-16TB75
		MT-16TB194	Anti-CD Toxin B mAb (clone TB194) (x1mg) / Paired with MT-16TB75
		MT-25TXB	C. difficile Toxin B recombinant protein (x1mg) (fragment without toxic activity)

8. For Rotavirus detection

		MT-16R15	Anti-Rotavirus mAb (clone R15) (x1mg) / Self Pairing!
		MT-25VP6	Rotavirus VP6 recombinant protein (x1mg)



1.2 Infectious diseases / Gastrointestinal

9. For Salmonella detection

9.1 For Salmonella paratyphi A detection

		MT-18SPA27	Anti-Salmonella paratyphi mAb (clone SPA27) (x1mg) / Self Pairing!
		MT-28SPAU	Inactivated Salmonella paratyphi A antigen (native extract) (x1mL)

9.2 For Salmonella typhi detection

		MT-18ST25	Anti-Salmonella typhi mAb (clone ST25) (x1mg) / Self Pairing!
		MT-28STU	Inactivated Salmonella typhi antigen (native extract) (x1mL)

9.3 For Salmonella paratyphi B and typhimurium detection

		MT-18SB2	Anti-Salmonella group B mAb (clone SB2) (x1mg) / Self Pairing!
		MT-28STMU	Inactivated Salmonella typhimurium antigen (native extract) (x1mL)
		MT-28SPBU	Inactivated Salmonella paratyphi B antigen (native extract) (x1mL)

9.4 For Salmonella enteritidis detection

		MT-18SE24	Anti-Salmonella enteritidis mAb (clone SE24) (x1mg) / Self Pairing!
		MT-28SEU	Inactivated Salmonella enteritidis antigen (native extract) (x1mL)

10. For Listeria detection

		MT-28LMU	Inactivated Listeria monocytogenes antigen (native extract) (x1mL)
--	---	----------	--

11. For Enterovirus detection

		MT-20EV	Anti-Enterovirus pAb (x1mg) / Self Pairing!
		MT-25ETV18	Enterovirus Echo18 recombinant protein (x1mg)
		MT-25ETV70	Enterovirus EV70 recombinant protein (x1mg)
		MT-25ETV71	Enterovirus EV71 recombinant protein (x1mg)
		MT-25ETVA16	Enterovirus CoxA16 recombinant protein (x1mg)
		MT-25ETVA24	Enterovirus CoxA24 recombinant protein (x1mg)
		MT-25ETVB3	Enterovirus CoxB3 recombinant protein (x1mg)

1.2 Infectious diseases / Gastrointestinal

12. For Entamoeba histolytica/dispar detection

		MT-16EH01	Anti-Entamoeba mAb (clone EH01) (x1mg) / Paired with MT-16EH11
		MT-16EH11	Anti-Entamoeba mAb (clone EH11) (x1mg) / Paired with MT-16EH01
		MT-25EHP	Entamoeba histolytica recombinant protein (x1mg)
		MT-25EDP	Entamoeba dispar recombinant protein (x1mg)

13. For Yersinia enterocolitica detection

		MT-28YE3U	Inactivated Yersinia enterocolitica O:3 antigen (native extract) (x1mL)
		MT-28YE9U	Inactivated Yersinia enterocolitica O:9 antigen (native extract) (x1mL)

14. For Shigella detection

		MT-28SFU	Inactivated Shigella flexneri antigen (native extract) (x1mL)
		MT-28SSU	Inactivated Shigella sonnei antigen (native extract) (x1mL)
		MT-28SDU	Inactivated Shigella dysenteriae antigen (native extract) (x1mL)
		MT-28SBU	Inactivated Shigella boydii antigen (native extract) (x1mL)

15. For E. coli O157 detection

15.1 For E. coli O157 detection

		MT-18E18	Anti-E. coli O157 mAb (clone E18) (x1mg) / Paired with MT-18E28
		MT-18E28	Anti-E. coli O157 mAb (clone E28) (x1mg) / Paired with MT-18E18
		MT-28EC7U	Inactivated E. coli O157 antigen (native extract) (x1mL)

115.2 For E. coli O157 Verotoxin 1 detection

		MT-25STX	E. coli O157 VT1 recombinant protein (x1mg)
--	---	----------	---

15.3 For E. coli O157 Verotoxin 2 detection

		MT-25VT2	E. coli O157 VT2 recombinant protein (x1mg)
--	---	----------	---

2. Antimicrobial Resistance

1. MCR-1

 MT-25MCR1	MCR1 recombinant protein (x1 mg)
---	----------------------------------

Carbapenemases

2. OXA-48

  MT-16OX21	Anti-OXA48 mAb (clone 21) (x1mg) / Paired with MT-16OX23
  MT-16OX23	Anti-OXA48 mAb (clone 23) (x1mg) / Paired with MT-16OX21
 MT-25OXA	OXA-48 recombinant protein (x1mg)

3. KPC

  MT-16KP05	Anti-KPC mAb (clone 05) (x1mg) / Paired with MT-16KP43
  MT-16KP43	Anti-KPC mAb (clone 43) (x1mg) / Paired with MT-16KP05 & MT-16KP58
  MT-16KP58	Anti-KPC mAb (clone 58) (x1mg) / Paired with MT-16KP43
 MT-25KPC	KPC recombinant protein (x1mg)

4. NDM

  MT-16ND10	Anti-NDM mAb (clone 10) (x1mg) / Paired with MT-16ND71
  MT-16ND22	Anti-NDM mAb (clone 22) (x1mg) / Paired with MT-16ND24
  MT-16ND24	Anti-NDM mAb (clone 24) (x1mg) / Paired with MT-16ND22
  MT-16ND71	Anti-NDM mAb (clone 71) (x1mg) / Paired with MT-16ND10
 MT-25NDM	NDM recombinant protein (x1mg)

5. VIM

  MT-16VI20	Anti-VIM mAb (clone 20) (x1mg) / Paired with MT-16VI21
  MT-16VI21	Anti-VIM mAb (clone 21) (x1mg) / Paired with MT-16VI20 & MT-16VI32
  MT-16VI32	Anti-VIM mAb (clone 32) (x1mg) / Paired with MT-16VI21
 MT-25VIM	VIM recombinant protein (x1mg)



_Monoclonal Antibodies



_Recombinant Proteins



_Paired with...

6. IMP

 	MT-16IM01	Anti-IMP mAb (clone 01) (x1mg) / Paired with MT-16IM05
 	MT-16IM05	Anti-IMP mAb (clone 05) (x1mg)/ Paired with MT-16IM01 & MT-16IM18
 	MT-16IM18	Anti-IMP mAb (clone 18) (x1mg)/ Paired with MT-16IM05
	MT-25IMP	IMP-1 recombinant protein (x1 mg)
	MT-25IMP4	IMP-4 recombinant protein (x1 mg)
	MT-25IMP8	IMP-8 recombinant protein (x1 mg)
	MT-25IMP13	IMP-13 recombinant protein (x1 mg)
	MT-25IMP14	IMP-14 recombinant protein (x1 mg)
	MT-25IMP22	IMP-22 recombinant protein (x1 mg)

3. Tumor & Inflammation markers

1. For human Calprotectin detection

 	MT-16CP14	Anti-Calprotectin mAb (clone CP14) (x1mg) / Self Pairing!
 	MT-16CP39	Anti-Calprotectin mAb (clone CP39) (x1mg) / Self Pairing!
 	MT-16CP43	Anti-Calprotectin mAb (clone CP43) (x1mg) / Paired with MT-16CP71
 	MT-16CP59	Anti-Calprotectin mAb (clone CP59) (x1mg) / Paired with MT-16CP71
  	MT-16CP71	Anti-Calprotectin mAb (clone CP71) (x1mg) / Self Pairing! / Paired with MT-16CP43 & MT-16CP59
	MT-25HCP	Human Calprotectin recombinant protein (x1mg)
	MT-25HCT	Human Calprotectin recombinant protein (x1mg)



Monoclonal Antibodies



Polyclonal Antibodies



Recombinant Proteins



Native Antigens Purified



Paired with...



Self Pairing

3. Tumor & Inflammation markers

2. For human Haemoglobin (FOB) detection

 	MT-16F22	Anti-Haemoglobin mAb (clone F22) (x1mg) / Paired with MT-16F52
 	MT-16F52	Anti-Haemoglobin mAb (clone F52) (x1mg) / Paired with MT-16F22
 	MT-20FOB	Anti-Haemoglobin pAb (x1mg) / Self Pairing!
	MT-29HHB	Human Haemoglobin protein (native extract) (x1mg)

3. For human Pancreatic Elastase detection

 	MT-16EL01	Anti-Elastase mAb (clone EL01) (x1mg) / Paired with MT-16EL04, MT-16EL40 & MT-16EL41
 	MT-16EL04	Anti-Elastase mAb (clone EL04) (x1mg) / Paired with MT-16EL01
 	MT-16EL40	Anti-Elastase mAb (clone EL40) (x1mg) / Paired with MT-16EL01
 	MT-16EL41	Anti-Elastase mAb (clone EL41) (x1mg) / Paired with MT-16EL01
 	MT-20HEM	Anti-Elastase pAb (x1mg) / Self Pairing!
	MT-25CELA3A	Human Elastase CELA3A recombinant protein (x1mg)
	MT-25CELA3B	Human Elastase CELA3B recombinant protein (x1mg)

4. Alfa-1-antitrypsin recombinant protein

	MT-25A1A	Alfa-1-antitrypsin recombinant protein (x1mg)
---	----------	---

5. For human Transferrin detection

 	MT-16TF8	Anti-Transferrin mAb (clone TF8) (x1mg) / Paired with MT-16TF16
 	MT-16TF16	Anti-Transferrin mAb (clone TF16) (x1mg) / Paired with MT-16TF8
	MT-29HTF	Human Transferrin protein (native extract) (x1mg)

6. For human IL-6 detection

	MT-25IL6	Human IL-6 recombinant protein (x1mg)
---	----------	---------------------------------------



Monoclonal Antibodies



Polyclonal Antibodies



Recombinant Proteins



Native Antigens Purified



Paired with...



Self Pairing

7. For human Lactoferrin detection

		MT-16LC20	Anti-Lactoferrin mAb (clone LC20) (x1mg) / Paired with MT-16LC37
		MT-16LC37	Anti-Lactoferrin mAb (clone LC37) (x1mg) / Paired with MT-16LC20
		MT-29HLF	Human Lactoferrin protein (native extract) (x1mg)

8. For human S100b detection

		MT-20HSBA	Anti-S100b pAb (affinity purified) (x1mg) / Self Pairing!
		MT-20HSBM	Anti-S100b pAb (protein A purified) (x1mg) / Self Pairing!
		MT-25HSB	Human S100b recombinant protein (x1mg)

9. For human S100A12 detection

		MT-16HS19	Anti-S100A12 mAb (clone HS19) (x1mg) / Paired with MT-16HS46
		MT-16HS46	Anti-S100A12 mAb (clone HS46) (x1mg) / Paired with MT-16HS19
		MT-20HSA	Anti-S100A12 pAb (affinity purified) (x1mg) / Self Pairing!
		MT-20HSM	Anti-S100A12 pAb (protein A purified) (x1mg) / Self Pairing!
		MT-25HSP	Human S100A12 recombinant protein (x1mg)

10. For human S100A9 and S100A8 detection

		MT-25HSA9	Human S100A9 recombinant protein (x1mg)
		MT-25HSA8	Human S100A8 recombinant protein (x1mg)

4. Other markers

1. Ferritin

		MT-25FER	Ferritin recombinant protein (x1 mg)
--	---	----------	--------------------------------------

Specificity | Sensitivity | Stability



Know How.

In vitro production



Storage:

- Frozen (-20°C)
- Refrigerated (2 - 8°C)



Deliveries.

Global service



Specific monoclonal screening:

- Fast
- Powerful
- Specific

“ Our expertise is to develop
accurate and reliable products
to improve your immunodiagnostic assays.

certest

RAW MATERIALS

Certest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, N°1
50840, San Mateo de Gállego, Zaragoza (Spain)
Tel. (+34) 976 520 354
RawMaterials@certest.es · www.certest.es

More *than*
diagnostics

Certest RawMaterials/GEN-1222EN

Modification rights reserved. All rights reserved. © Certest Biotec, S.L.
The products, services and data set out in this document may suffer changes and/or variations on the texts and pictures shown.