

VIASURE

by **certest**

# VIASURE Universal Pathogens Extraction Kit for VIASURE V-Quick32

(VS-EAU0196VQ)

**IVDR**  
CERTIFIED



**Ready to use reagents**  
with preloaded plates.



**Stable at Room  
Temperature.**



**Magnetic beads**  
technology.



**Purification of viral,  
bacterial, parasitic and  
fungal DNA and RNA.**



**Compatible with a wide  
range of sample types:**

- Swab specimens  
(i.e. respiratory tract,  
urogenital tract and  
lesions)
- Saliva
- Sputum
- BAL & BAS
- Feces
- Urine
- Semen
- Whole blood
- Blood culture
- Plasma
- Serum
- Cerebrospinal fluid (CSF)
- Gastric biopsies
- Isolated culture colony



VIASURE/EK-UNIVVQUICK-EN

Modification rights reserved. All rights reserved. © Certest Biotech, S.L.

The products, services and data set out in this document may suffer changes and/or variations on the texts and pictures shown.

**certest**

Certest Biotech, S.L.  
Pol. Industrial Río Gállego II - Calle J, N°1 50840, San  
Mateo de Gállego, Zaragoza (Spain)  
Tel. (+34) 976 520 354 | viasure@certest.es  
[www.certest.es](http://www.certest.es)