

Technical Sheet

Product Name:	Helicobacter pylori protein
Business Unit:	CerTest Raw Materials

TECHNICAL DETAILS:

Code:	MT-24PYL
Physical State:	Liquid (frozen)
Source:	<i>H. pylori whole cell sonicate</i>
Description:	H. pylori protein
Storage Conditions:	Shipping: -20°C. Storage at -20°C. Avoid multiple/freeze thaw cycles by storing multiple aliquots at -20°C.
Concentration	1 mg/ml
Health & Safety Information:	This material has been centrifuged to remove whole bacteria and other insoluble material. It is a soluble protein containing urease, cagA & cytotoxin A antigens. This product contains sodium azide as preservative. Although the amount is 0.09%, this substance is poisonous and hazardous, and it should be handled by trained staff only. This material has been centrifuged to remove whole bacteria and other insoluble material. Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.
Presentation:	PBS (1x) buffer with additives and salts, pH 7.4, 0.09% Sodium Azide (NaN ₃).

Date of Technical sheet: **22/05/2024**

Signed:



Rafael Clavería
Bioscience Manager

Code of the document and version: **Tshb-MT-24PYL v.00**

CHANGE CONTROL:

CHANGE CONTROL		
Version nº	Changes	Date
00	Original version of the file	22/05/2024