


# Rapid Test

One-step immunochromatographic tests

Fast | Easy | Reliable | Simple



Certest Immunodiagnostic is committed to providing a comprehensive portfolio of diagnostic reagents for detecting infectious diseases and different markers, covering different needs of formats, types of use, throughput, precision, and time from sample to result.


### IMMUNOCHROMATOGRAPHY

provides rapid test solutions for aiding the professional of diagnostic with an immediate solution, needing no equipment, and having results usually in less than 15 minutes.

### TURBIDIMETRY

provides diagnostics solutions for automated and quantitative results and higher throughputs.

**CLIA** provides a high sensitivity, automated diagnostic solution based on the chemiluminescent amplification of an immunocapture.



## Rapid Test

Our products are the result of years of experience detecting the presence of antigens with the lateral flow **immunochromatographic** technique.

The devices are based on the coloured latex particles. Showing results on bands of different colours facilitates the interpretation of the result.


**Certest Lateral Flow Immunochromatography** portfolio mainly covers the diagnostic of the gastrointestinal, respiratory infectious diseases and the detection of inflammatory or tumor markers in the Gastrointestinal tract.

- ✓ **Single device**, multiple tests from one sample.
- ✓ **One sample**, one preparation, panel testing.
- ✓ Multiple analytes diagnosed in **10 minutes**.
- ✓ **Reduced risk** of sample and result mishandling.
- ✓ Operational efficiency.



### Multianalyte COMBO (1, 2, 3 or 4 strips)

*Our cut-off tests are equipped with special vials which facilitates fecal sample collection.*



| Fecal antigens                                                                      |                                                                                                                                                                                                                                                                                                         | Respiratory antigens                                                                                                                                                       |  |  |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|
|    |  <b>VIRUSES</b><br>Adenovirus<br>Astrovirus<br>Enterovirus<br>Norovirus GI/GII<br>Rotavirus                                                                                                                            |  <b>VIRUSES</b><br>Flu A<br>RSV<br>Resp. Adenovirus<br>SARS-CoV-2                         |  |  |
|                                                                                     |  <b>BACTERIA</b><br>Helicobacter pylori<br>Clostridium difficile antigen GDH<br>Escherichia coli O157<br>Campylobacter<br>Salmonella<br>Salmonella Typhi<br>Salmonella Paratyphi<br>Yersinia O:3<br>Yersinia O:9       |  <b>BACTERIA</b><br>Strep. A                                                              |  |  |
|                                                                                     |  <b>PARASITES</b><br>Cryptosporidium parvum<br>Giardia lamblia<br>Entamoeba                                                                                                                                            |                                                                                                                                                                            |  |  |
|                                                                                     |                                                                                                                                                                                                                                                                                                         |                                                                                                                                                                            |  |  |
|  |  <b>VIRUSES</b><br>Rotavirus + Adenovirus<br>Norovirus GI + Norovirus GII<br>Rotavirus + Norovirus                                                                                                                   |  <b>VIRUSES</b><br>Flu A + Flu B<br>RSV + Resp. Adenovirus<br>SARS-CoV-2 + Flu A        |  |  |
|                                                                                     |  <b>BACTERIA</b><br>H. pylori + Transferrin<br>Clostridium difficile (GDH) + Lactoferrin<br>Clostridium difficile Toxin A + B<br>EHEC VT1 + VT2<br>Salmonella typhi + paratyphi<br>Yersinia enterocolitica O:3 + O:9 |                                                                                                                                                                            |  |  |
|                                                                                     |  <b>PARASITES</b><br>Crypto + Giardia                                                                                                                                                                                |                                                                                                                                                                            |  |  |
|                                                                                     |                                                                                                                                                                                                                                                                                                         |                                                                                                                                                                            |  |  |
|  |  <b>VIRUSES</b><br>Rotavirus + Adenovirus + Astrovirus<br>Rotavirus + Adenovirus + Norovirus                                                                                                                         |  <b>VIRUSES</b><br>Flu A + Flu B + RSV<br>SARS-CoV-2 + Flu A + Flu B                    |  |  |
|                                                                                     |  <b>BACTERIA</b><br>Clostridium difficile GDH + Toxin A + Toxin B<br>E. coli O157 + VT1 + VT2                                                                                                                        |                                                                                                                                                                            |  |  |
|                                                                                     |  <b>PARASITES</b><br>Crypto + Giardia + Entamoeba                                                                                                                                                                    |                                                                                                                                                                            |  |  |
|                                                                                     |                                                                                                                                                                                                                                                                                                         |                                                                                                                                                                            |  |  |
|  |  <b>VIRUSES</b><br>Rotavirus + Adenovirus + Astrovirus + Norovirus<br>GDH + Toxin A + Toxin B + Lactoferrin                                                                                                          |  <b>VIRUSES</b><br>Flu A + Flu B + RSV + Adenovirus<br>SARS-CoV-2 + Flu A + Flu B + RSV |  |  |
|                                                                                     |  <b>BACTERIA</b><br>C. Difficile GDH + Toxin A + Toxin B + Lactoferrin                                                                                                                                               |                                                                                                                                                                            |  |  |

**Inflammatory and  
tumor markers**



 **MARKERS**

Hemoglobin (FOB)  
Transferrin  
Calprotectin  
Lactoferrin

**Urine**

 **BACTERIA**

Legionella  
Streptococcus pneumoniae



 **MARKERS**

FOB 50 + 200  
FOB + Transferrin  
H. pylori + Transferrin  
Calprotectin 50 + 200  
Calprotectin + Lactoferrin  
Clostridium difficile (GDH) + Lactoferrin

 **BACTERIA**

Streptococcus pneumoniae + Legionella



 **MARKERS**

FOB + Transferrin + Calprotectin + Lactoferrin  
GDH + Toxin A + Toxin B + Lactoferrin

**Our tests have two  
differentiated lines.**

Control line - green color  
Positive line - red color

**Possibility to add *positive*  
and *negative* controls.**





Fecal Antigens

 GASTROINTESTINAL VIRUS



ADENOVIRUS

| Reference  | Product name                                                       |
|------------|--------------------------------------------------------------------|
| AD820001V  | Adenovirus (20 test cassettes with stool vials)                    |
| AD820001VC | Adenovirus (20 test cassettes with stool vials + Positive Control) |
| AD220001P  | Adenovirus 20 Positive Controls Kit                                |



ASTROVIRUS

| Reference  | Product name                                                       |
|------------|--------------------------------------------------------------------|
| AT820001V  | Astrovirus (20 test cassettes with stool vials)                    |
| AT820001VC | Astrovirus (20 test cassettes with stool vials + Positive Control) |
| AT220001P  | Astrovirus 20 Positive Controls Kit                                |



ENTEROVIRUS

| Reference  | Product name                                                        |
|------------|---------------------------------------------------------------------|
| EV820001V  | Enterovirus (20 test cassettes with stool vials)                    |
| EV820001VC | Enterovirus (20 test cassettes with stool vials + Positive Control) |
| EV220001P  | Enterovirus 20 Positive Controls Kit                                |



NOROVIRUS

| Reference  | Product name                                                      |
|------------|-------------------------------------------------------------------|
| NN820001V  | Norovirus (20 test cassettes with stool vials)                    |
| NN820001VC | Norovirus (20 test cassettes with stool vials + Positive Control) |
| NN220001P  | Norovirus 20 Positive Controls Kit                                |



ROTAVIRUS

| Reference | Product name                                                             |
|-----------|--------------------------------------------------------------------------|
| R820001V  | <b>Rotavirus</b> (20 test cassettes with stool vials)                    |
| R820001VC | <b>Rotavirus</b> (20 test cassettes with stool vials + Positive Control) |
| R220001P  | <b>Rotavirus</b> 20 Positive Controls Kit                                |



COMBO

| Reference  | Product name                                                                                  |
|------------|-----------------------------------------------------------------------------------------------|
| XR882001V  | <b>Rotavirus + Adenovirus</b> (20 test combo cassettes with stool vials)                      |
| XR882001VC | <b>Rotavirus + Adenovirus</b> (20 test combo cassettes with stool vials + Positive Controls)  |
| NV882001V  | <b>Norovirus GI + GII</b> (20 test combo cassettes with stool vials)                          |
| NV882001VC | <b>Norovirus GI + GII</b> (20 test combo cassettes with stool vials + Positive Controls)      |
| RN882001V  | <b>Rotavirus + Noro GI/GII</b> (20 test combo cassettes with stool vials)                     |
| RN882001VC | <b>Rotavirus + Noro GI/GII</b> (20 test combo cassettes with stool vials + Positive Controls) |



TRIPLE

| Reference  | Product name                                                                                               |
|------------|------------------------------------------------------------------------------------------------------------|
| RT872001V  | <b>Rotavirus + Adenovirus + Astrovirus</b> (20 test triple cassettes with stool vials)                     |
| RT872001VC | <b>Rotavirus + Adenovirus + Astrovirus</b> (20 test triple cassettes with stool vials + Positive Controls) |
| RA872001V  | <b>Rotavirus + Adenovirus + Norovirus</b> (20 test triple cassettes with stool vials)                      |
| RA872001VC | <b>Rotavirus + Adenovirus + Norovirus</b> (20 test triple cassettes with stool vials + Positive Controls)  |



QUAD

| Reference  | Product name                                                                                     |
|------------|--------------------------------------------------------------------------------------------------|
| RG862001V  | <b>Rota + Adeno + Astro + Noro</b> (20 test quad cassettes with stool vials)                     |
| RG862001VC | <b>Rota + Adeno + Astro + Noro</b> (20 test quad cassettes with stool vials + Positive Controls) |



## Fecal Antigens

### GASTROINTESTINAL BACTERIA



#### HELICOBACTER PYLORI

| Reference | Product name                                                             |
|-----------|--------------------------------------------------------------------------|
| P820001V  | <b>H. Pylori</b> (20 test cassettes with stool vials)                    |
| P820001VC | <b>H. Pylori</b> (20 test cassettes with stool vials + Positive Control) |
| P220001P  | <b>H. Pylori</b> 20 Positive Controls Kit                                |



#### CLOSTRIDIUM DIFFICILE

| Reference  | Product name                                                                                     |
|------------|--------------------------------------------------------------------------------------------------|
| GD820001V  | <b>Clostridium Difficile antigen GDH</b> (20 test cassettes with stool vials)                    |
| GD820001VC | <b>Clostridium Difficile antigen GDH</b> (20 test cassettes with stool vials + Positive Control) |
| GD220001P  | <b>Clostridium Difficile antigen GDH</b> 20 Positive Controls Kit                                |
| TA220001P  | <b>Clostridium Difficile Toxin A</b> 20 Positive Controls Kit                                    |
| TB220001P  | <b>Clostridium Difficile Toxin B</b> 20 Positive Controls Kit                                    |



#### ESCHERICHIA COLI

| Reference | Product name                                                                   |
|-----------|--------------------------------------------------------------------------------|
| E820001V  | <b>E. Coli O157:H7</b> (20 test cassettes with stool vials)                    |
| E820001VC | <b>E. Coli O157:H7</b> (20 test cassettes with stool vials + Positive Control) |
| E220001P  | <b>E. Coli O157:H7</b> 20 Positive Controls Kit                                |
| VT220001P | <b>E. Coli O157:H7</b> Verotoxin 1, 20 Positive Controls Kit                   |
| VX220001P | <b>E. Coli O157:H7</b> Verotoxin 2, 20 Positive Controls Kit                   |



#### CAMPYLOBACTER

| Reference  | Product name                                                                 |
|------------|------------------------------------------------------------------------------|
| CA820001V  | <b>Campylobacter</b> (20 test cassettes with stool vials)                    |
| CA820001VC | <b>Campylobacter</b> (20 test cassettes with stool vials + Positive Control) |
| CA220001P  | <b>Campylobacter</b> 20 Positive Controls Kit                                |

SALMONELLA

| Reference  | Product name                                                              |
|------------|---------------------------------------------------------------------------|
| SA820001V  | <b>Salmonella</b> (20 test cassettes with stool vials)                    |
| SA820001VC | <b>Salmonella</b> (20 test cassettes with stool vials + Positive Control) |
| SA220001P  | <b>Salmonella</b> 20 Positive Controls Kit                                |

SALMONELLA TYPHI

| Reference  | Product name                                                                    |
|------------|---------------------------------------------------------------------------------|
| ST820001V  | <b>Salmonella Typhi</b> (20 test cassettes with stool vials)                    |
| ST820001VC | <b>Salmonella Typhi</b> (20 test cassettes with stool vials + Positive Control) |
| ST220001P  | <b>Salmonella Typhi</b> 20 Positive Controls Kit                                |

SALMONELLA PARATYPHI

| Reference  | Product name                                                                        |
|------------|-------------------------------------------------------------------------------------|
| SP820001V  | <b>Salmonella Paratyphi</b> (20 test cassettes with stool vials)                    |
| SP820001VC | <b>Salmonella Paratyphi</b> (20 test cassettes with stool vials + Positive Control) |
| SP220001P  | <b>Salmonella Paratyphi</b> 20 Positive Controls Kit                                |

YERSINISA O:3

| Reference  | Product name                                                                               |
|------------|--------------------------------------------------------------------------------------------|
| YE820001V  | <b>Yersinia Enterocolitica O:3</b> (20 test cassettes with stool vials)                    |
| YE820001VC | <b>Yersinia Enterocolitica O:3</b> (20 test cassettes with stool vials + Positive Control) |
| YE220001P  | <b>Yersinia Enterocolitica O:3</b> 20 Positive Controls Kit                                |

YERSINISA O:9

| Reference  | Product name                                                                               |
|------------|--------------------------------------------------------------------------------------------|
| YC820001V  | <b>Yersinia Enterocolitica O:9</b> (20 test cassettes with stool vials)                    |
| YC820001VC | <b>Yersinia Enterocolitica O:9</b> (20 test cassettes with stool vials + Positive Control) |
| YC220001P  | <b>Yersinia Enterocolitica O:9</b> 20 Positive Controls Kit                                |



Fecal Antigens



COMBO

| Reference   | Product name                                                                                                  |
|-------------|---------------------------------------------------------------------------------------------------------------|
| CD882001V   | <b>Clostridium Difficile Toxin A + Toxin B</b> (20 test combo cassettes with stool vials)                     |
| CD882001VC  | <b>Clostridium Difficile Toxin A + Toxin B</b> (20 test combo cassettes with stool vials + Positive Controls) |
| PT882001VF  | <b>H. Pylori + Transferrin</b> (20 test combo cassettes with stool vials)                                     |
| PT882001VFC | <b>H. Pylori + Transferrin</b> (20 test combo cassettes with stool vials + Positive Controls)                 |
| GF882001VF  | <b>GDH + Lactoferrin</b> (20 test combo cassettes with stool vials)                                           |
| GF882001VFC | <b>GDH + Lactoferrin</b> (20 test combo cassettes with stool vials + Positive Controls)                       |
| SX882001V   | <b>EHEC VT1 + VT2</b> (20 test combo cassettes with stool vials)                                              |
| SX882001VC  | <b>EHEC VT1 + VT2</b> (20 test combo cassettes with stool vials + Positive Controls)                          |
| TP882001V   | <b>Salmonella Typhi + Paratyphi</b> (20 test combo cassettes with stool vials)                                |
| TP882001VC  | <b>Salmonella Typhi + Paratyphi</b> (20 test combo cassettes with stool vials + Positive Controls)            |
| YS882001V   | <b>Yersinia Enterocolitica O:3 + O:9</b> (20 test combo cassettes with stool vials)                           |
| YS882001VC  | <b>Yersinia Enterocolitica O:3 + O:9</b> (20 test combo cassettes with stool vials + Positive Controls)       |



TRIPLE

| Reference  | Product name                                                                                                |
|------------|-------------------------------------------------------------------------------------------------------------|
| GX872001V  | <b>C. Difficile GDH + Toxin A + Toxin B</b> (20 test triple cassettes with stool vials)                     |
| GX872001VC | <b>C. Difficile GDH + Toxin A + Toxin B</b> (20 test triple cassettes with stool vials + Positive Controls) |
| EX872001V  | <b>E. Coli O157:H7 + VT1 + VT2</b> (20 test triple cassettes with stool vials)                              |
| EX872001VC | <b>E. Coli O157:H7 + VT1 + VT2</b> (20 test triple cassettes with stool vials + Positive Controls)          |



QUAD

| Reference   | Product name                                                                                                            |
|-------------|-------------------------------------------------------------------------------------------------------------------------|
| LX862001VF  | <b>C. Difficile GDH + Toxin A + Toxin B + Lactoferrin</b> (20 test quad cassettes with stool vials)                     |
| LX862001VFC | <b>C. Difficile GDH + Toxin A + Toxin B + Lactoferrin</b> (20 test quad cassettes with stool vials + Positive Controls) |

 **GASTROINTESTINAL PARASITES**

 **CRYPTOPROSPORIDIUM**

| Reference | Product name                                                          |
|-----------|-----------------------------------------------------------------------|
| K820001V  | <b>Crypto</b> (20 test cassettes with stool vials)                    |
| K820001VC | <b>Crypto</b> (20 test cassettes with stool vials + Positive Control) |
| K220001P  | <b>Crypto</b> 20 Positive Controls Kit                                |

 **GIARDIA**

| Reference  | Product name                                                            |
|------------|-------------------------------------------------------------------------|
| GW820001V  | <b>Giardia+</b> (20 test cassettes with stool vials)                    |
| GW820001VC | <b>Giardia+</b> (20 test cassettes with stool vials + Positive Control) |
| GW220001P  | <b>Giardia+</b> 20 Positive Controls Kit                                |

 **ENTAMOEBA**

| Reference  | Product name                                                             |
|------------|--------------------------------------------------------------------------|
| EH820001V  | <b>Entamoeba</b> (20 test cassettes with stool vials)                    |
| EH820001VC | <b>Entamoeba</b> (20 test cassettes with stool vials + Positive Control) |
| EH220001P  | <b>Entamoeba</b> 20 Positive Controls Kit                                |

 **COMBO**

| Reference  | Product name                                                                           |
|------------|----------------------------------------------------------------------------------------|
| NW882001V  | <b>Crypto + Giardia</b> (20 test combo cassettes with stool vials)                     |
| NW882001VC | <b>Crypto + Giardia</b> (20 test combo cassettes with stool vials + Positive Controls) |

 **TRIPLE**

| Reference  | Product name                                                                                |
|------------|---------------------------------------------------------------------------------------------|
| NE872001V  | Crypto + Giardia + Entamoeba (20 test combo cassettes with stool vials)                     |
| NE872001VC | Crypto + Giardia + Entamoeba (20 test combo cassettes with stool vials + Positive Controls) |



## Respiratory Antigens



### RESPIRATORY VIRUS



#### INFLUENZA

| Reference | Product name                                                        |
|-----------|---------------------------------------------------------------------|
| Y820001PC | Influenza A (20 test cassettes with stool vials + Positive Control) |
| Y220001P  | Influenza A 20 Positive Controls Kit                                |
| YB220001P | Influenza B 20 Positive Controls Kit                                |



#### RSV (RESPIRATORY SYNCYTIAL VIRUS)

| Reference  | Product name                                                |
|------------|-------------------------------------------------------------|
| RV820001PC | RSV (20 test cassettes with stool vials + Positive Control) |
| RV220001P  | RSV 20 Positive Controls Kit                                |



#### ADENOVIRUS RESPIRATORY

| Reference  | Product name                                                             |
|------------|--------------------------------------------------------------------------|
| AP820001PC | Adenovirus Resp. (20 test cassettes with stool vials + Positive Control) |
| AP220001P  | Adenovirus Resp. 20 Positive Controls Kit                                |



#### SARS COV-2

| Reference  | Product name                                                       |
|------------|--------------------------------------------------------------------|
| SC820001P  | SARS CoV-2 (20 test cassettes with stool vials)                    |
| SC820001PC | SARS CoV-2 (20 test cassettes with stool vials + Positive Control) |
| SC860001PC | SARS CoV-2 (60 test cassettes with stool vials + Positive Control) |
| SC220001P  | SARS CoV-2 20 Positive Controls Kit                                |



## COMBO

| Reference  | Product name                                                                          |
|------------|---------------------------------------------------------------------------------------|
| Z882001PC  | <b>Influenza A + B</b> (20 test cassettes with stool vials + Positive Control)        |
| VR882001PC | <b>RSV + Adenovirus Resp.</b> (20 test cassettes with stool vials + Positive Control) |
| SZ882001PC | <b>SARS CoV-2 + Flu A</b> (20 test cassettes with stool vials + Positive Control)     |



## TRIPLE

| Reference  | Product name                                                                              |
|------------|-------------------------------------------------------------------------------------------|
| ZS872001PC | <b>Influenza A + B + RSV</b> (20 test cassettes with stool vials + Positive Control)      |
| SS872001PC | <b>SARS CoV-2 + Flu A + Flu B</b> (20 test cassettes with stool vials + Positive Control) |



## QUAD

| Reference  | Product name                                                                                       |
|------------|----------------------------------------------------------------------------------------------------|
| ZV862001PC | <b>Influenza A + B + RSV + Adeno Resp.</b> (20 test cassettes with stool vials + Positive Control) |
| SV862001PC | <b>SARS CoV-2 + Flu A + Flu B + RSV</b> (20 test cassettes with stool vials + Positive Control)    |



## RESPIRATORY BACTERIA

### STREP. A

| Reference | Product name                                                                                          |
|-----------|-------------------------------------------------------------------------------------------------------|
| T820201PC | <b>Strep A</b> (20 test cassettes with diluent bottles + tubes and pipettes vials + Positive Control) |
| T220001P  | <b>Strep A</b> 20 Positive Controls Kit                                                               |



Any test KIT which includes positive control could have negative control too adding C at the end of the reference.

| Reference | Product name                              |
|-----------|-------------------------------------------|
| CN320001P | Negative Control 20 Negative Controls Kit |



Urine Antigens



URINE BACTERIA



LEGIONELLA PNEUMOPHILA

| Reference  | Product name                                                                                    |
|------------|-------------------------------------------------------------------------------------------------|
| LP820201PC | <b>Legionella</b> (20 test cassettes with diluent bottle + Positive Control + Negative Control) |
| LP220001P  | <b>Legionella</b> 20 Positive Controls Kit                                                      |



STREPTOCOCCUS PNEUMONIAE

| Reference  | Product name                                                                                       |
|------------|----------------------------------------------------------------------------------------------------|
| SN820201PC | <b>S. pneumoniae</b> (20 test cassettes with diluent bottle + Positive Control + Negative Control) |
| SN220001P  | <b>S. pneumoniae</b> 20 Positive Controls Kit                                                      |



COMBO

| Reference  | Product name                                                                                          |
|------------|-------------------------------------------------------------------------------------------------------|
| NL882201PC | <b>S. pneumoniae + Legionella</b> (20 test cassettes with diluent bottle + Positive Controls + N. C.) |



Inflammatory & Tumor Makers



FOB

| Reference | Product name                                                       |
|-----------|--------------------------------------------------------------------|
| F820001F  | <b>FOB</b> (20 test cassettes with stool vials)                    |
| F820001FC | <b>FOB</b> (20 test cassettes with stool vials + Positive Control) |
| F220001P  | <b>FOB</b> 20 Positive Controls Kit                                |



CALPROTECTIN

| Reference  | Product name                                                                |
|------------|-----------------------------------------------------------------------------|
| CP820001F  | <b>Calprotectin</b> (20 test cassettes with stool vials)                    |
| CP820001FC | <b>Calprotectin</b> (20 test cassettes with stool vials + Positive Control) |
| CP220001P  | <b>Calprotectin</b> 20 Positive Controls Kit                                |



TRANSFERRIN

| Reference  | Product name                                                               |
|------------|----------------------------------------------------------------------------|
| TF820001F  | <b>Transferrin</b> (20 test cassettes with stool vials)                    |
| TF820001FC | <b>Transferrin</b> (20 test cassettes with stool vials + Positive Control) |
| TF220001P  | <b>Transferrin</b> 20 Positive Controls Kit                                |



LACTOFERRIN

| Reference  | Product name                                                               |
|------------|----------------------------------------------------------------------------|
| LT820001F  | <b>Lactoferrin</b> (20 test cassettes with stool vials)                    |
| LT820001FC | <b>Lactoferrin</b> (20 test cassettes with stool vials + Positive Control) |
| LT220001P  | <b>Lactoferrin</b> 20 Positive Controls Kit                                |



COMBO

| Reference  | Product name                                                                                    |
|------------|-------------------------------------------------------------------------------------------------|
| FB882001F  | <b>FOB 50 + 200</b> (20 test combo cassettes with stool vials)                                  |
| FB882001FC | <b>FOB 50 + 200</b> (20 test combo cassettes with stool vials + Positive Controls)              |
| CC882001F  | <b>Calprotectin 50 + 200</b> (20 test combo cassettes with stool vials)                         |
| CC882001FC | <b>Calprotectin 50 + 200</b> (20 test combo cassettes with stool vials + Positive Controls)     |
| FT882001F  | <b>FOB + Transferrin</b> (20 test combo cassettes with stool vials)                             |
| FT882001FC | <b>FOB + Transferrin</b> (20 test combo cassettes with stool vials + Positive Controls)         |
| CL882001F  | <b>Calprotectin + Lactoferin</b> (20 test combo cassettes with stool vials)                     |
| CL882001FC | <b>Calprotectin + Lactoferin</b> (20 test combo cassettes with stool vials + Positive Controls) |



QUAD

| Reference  | Product name                                                                                                         |
|------------|----------------------------------------------------------------------------------------------------------------------|
| FC862001F  | <b>FOB + Transferrin + Calprotectin + Lactoferrin</b> (20 test combo cassettes with stool vials)                     |
| FC862001FC | <b>FOB + Transferrin + Calprotectin + Lactoferrin</b> (20 test combo cassettes with stool vials + Positive Controls) |

## Advantages:



Noninvasive diagnostic



All included. No additional equipment needed



Low operational costs



Immediate results



Ease of use and interpretation

**Certest Immunodiagnostic** offers simple, sensitive, specific screening assays for a presumptive diagnosis of infectious diseases.



Find us on

@certestbiotec

**certest**

IMMUNODIAGNOSTIC

**Certest Biotec, S.L.**

Pol. Industrial Río Gállego II · Calle J, N°1

50840, San Mateo de Gállego, Zaragoza (Spain)

Tel. (+34) 976 520 354

[salesidx@certest.es](mailto:salesidx@certest.es)

[www.certest.es](http://www.certest.es)

Certest Rapid Test/GEN-0624EN

Modification rights reserved. All rights reserved. © Certest Biotec, S.L.  
The products, services and data set out in this document may suffer changes and/or variations on the texts and pictures shown.